



ISO 9001 certified



BIONOTE

Product Catalog ver3.0

CIA 03-03



Diagnostic Test Kits for Industrial Animals

Rapid

- [Avian]**
 - AIV Ag
 - H5 AIV Ag
 - AIV/H5 AIV Ag
 - AIV Ab
 - NDV Ag
 - IBDV Ag
- [Bovine]**
 - B.brucella Ab
 - Bovine TB Ab
 - FMD NSP Ab
 - Rota Ag
- [Swine]**
 - PED Ag
 - TGE Ag
 - TGE/PED Ag

ELISA

- [Avian]**
 - AIV Ab
 - H5 AIV Ab
 - NDV velo Ab
- [Bovine]**
 - B.brucella Ab
 - BTB Ab
 - FMD NSP Ab
- [Swine]**
 - CSFV Ab
 - PRRS Ab

PCR

- [Avian]**
 - General PCR Kit
 - AIV Detection Kit
 - NDV, IBDV Detection Kit
- [Real-Time PCR Kit]**
 - AIV Detection Kit
 - Other poultry disease virus Detection Kit
- [Swine]**
 - Real-Time PCR kit
 - PRRS, CSFV, PCV, Detection Kit
 - SIV Detection Kit





Global expertise in In Vitro Diagnostics

BIONOTE Product Catalog(Ver.3.0)

BIONOTE was established at the beginning of 2003, and is considered a expert of In –Vitro Diagnostics for veterinary in needs. We deliver inspirations, innovation and quality solutions for our customers and consumers to improve animal health care. BIONOTE has own automated facility to manufacture a wide range of high value products developed by high Qualified R&D center. BIONOTE manufactures diagnostic products to accordance with ISO9001:ISO13485 certification, and has expanded overseas sales network over 90 countries and is still continuing to grow.

Company History

- * 2009-2011 Expanding business**

2011.01 BIONOTE was merged into with Allere.

2009.03 Moved to current facility at 2-9,Seogu-dong, Hwaseong-si, Gyeonggi-do, Korea 445-170

2009.02 Company name was changed from Animal Genetics, Inc. to BioNote, Inc.
- * 2007-2008 Take – off stage**

2007.12 Won a prize at the 2007 self-audit evaluation(NVRQS)

2007.07 Certificated as a promising small & medium business of Gyeonggi province
- * 2005-2006 Installing Base profit**

2006.06 Certificated as the venture company

2006.01 Acquired new manufacturing facility area in Hwaseong-si(City)


2005.08 Moved to new manufacturing facility
- * 2003-2004 Improving Business Environment**

2004.07 Acquired technical innovation subject of small& medium business

2003.12 Acquired certificate of manufacture for animal health product [NVRQS]

2003.11 Acquired certificate of import for animal product [NVRQS]

2003.03 Establish "Animal Genetics, Inc."



Anigen Rapid AIV Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of avian influenza type A viruses antigen (H1~H16) in avian cloaca, trachea, kidney swab and feces.

Background

Influenza viruses that infect birds are called "Avian Influenza virus (AIV)" with only influenza type A viruses infecting birds. Influenza type A viruses are divided into subtypes based on two proteins, hemagglutinin (HA) and neuraminidase (NA), on the surface of the virus. AIV can infect chickens, turkeys, pheasants, quail, ducks, geese, and guinea fowl, as well as a wide variety of other birds. AIV is spread primarily through direct contact from infected birds to healthy birds, and through indirect contact with contaminated equipment and materials. The virus is excreted through infected birds' feces and secretions from their noses, mouths, and eyes.


Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - Monoclonal anti-AIV(Capture) - **AIV Type A** - Monoclonal anti-AIV(Detector)
 - Monoclonal antibody against Type A specific nucleoprotein of Avian influenza virus
- Purpose: Detection of Avian influenza type A virus(H1~H16)
- Specimen: Avian cloaca, trachea, kidney or feces swab
- Reading time: 5~10 minutes
- Sensitivity: 100% by farm (n=19), 77.3% by feces (n=150)
- Specificity: 100% vs. HA, PCR (n=1,402)
- No cross reaction against NDV, IBV, ILTV, Pneumovirus, Reovirus, IBDV, MDV, Mycoplasma and REV
- Detection limit: $10^{3.1}$ EID₅₀/ml
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 30 multi-devices/Kit, 10 and 25 individual devices/Kit
- The Pigeons have very low possibility of AIV infection and the specimens from pigeon feces may show unusual result in the kit.
To test the pigeon feces, Choose Anigen pigeon feces preparation kit that is commercially available.

Special Features

- Optimal screening method of AIV in the farm or laboratory
- Detection of all avian influenza type A viruses (H1N1 ~ H16N9)
- High Sensitivity compared with PCR, virus isolation method, HA and other immunoassays
- High specificity: No cross-reaction with other avian viruses.
- Detection of AIV in Chickens, Ducks, Turkeys, and Geese
- World's first commercialized rapid test kit for detection of avian influenza virus since December 2003
- Easy assay procedure: No technical expertise required, and Saving Labor & Time
- No additional equipment is required.
- It is recommended to test 10 birds in every AIV suspected flock for the excellent result.

Test Procedures



1. Collect a cloaca feces swab, tracheal swab, kidney swab or scattered wet feces swab.
2. Dilute the swabbed sample with assay diluent and drop 4 drops.
3. Wait for 5~10 minutes and read the test result.



CAT.No.: RG15-01

Performance Data

1. Evaluation of Influenza virus type A Reactivity

A/RP/8/34 (H1N1)	Positive
A/Singapore/1/57 (H2N2)	Positive
A/Duck/Ukraine/1/63 (H3N8)	Positive
A/Duck/Czechoslovakia/56 (H4N6)	Positive
A/Duck/Vietnam/TG24-01/05 (H5N1)	Positive
A/Duck/Hong Kong/820/80 (H5N3)	Positive*
A/Shearwater/Australia/1/72 (H6N5)	Positive
A/Duck/Hong Kong/301/78 (H7N1)	Positive
A/chicken/Germany/R28/03 (H7N7)	Positive*
A/Turkey/Ontario/6118/68 (H8N4)	Positive
A/Turkey/Wisconsin/1/66 (H9N2)	Positive
A/Chick/Germany/N/49 (H10N7)	Positive
A/Duck/England/56 (H11N6)	Positive
A/Duck/Alberta/60/76 (H12N5)	Positive
A/Gull/Maryland/707/77 (H13N6)	Positive
A/Mallard/Gurjev/263/82 (H14N5)	Positive
A/Shearwater/West Australia/2576/79 (H15N9)	Positive
A/BHG/Sweden/5/99 (H16N3)	Positive**

Evaluated at National Veterinary Research & Quarantine service (NVRQS), Korea
 *Evaluated at Friedrich Loeffler Institute (FLI), Germany (OIE AIV Ref. Lab)
 **Evaluated at CSIRO, Australian Animal Health Laboratory (OIE AIV Ref. Lab.)

2. Detection Limit Study

Specimen	AIV A/ch/Germ/R28/03 (H7N7)		
	14/3/R 5*103.8EID60/ml	14/1/R 5*103.1EID50/ml	14/2/K 5*102.8EID50/ml
Anigen AIV Ag	Strong positive	Strong positive	Negative

Evaluated at Friedrich Loeffler Institute (FLI), Germany (OIE AIV Ref. Lab)

3. Early Detection Study

Challenge method : SPF Chicken (H9N2, 106EID50/0.1ml) by nasal inoculation

	Detection Rate(%) by Day Post Inoculation						
	1 DPI	2 DPI	3 DPI	4 DPI	5 DPI	7 DPI	10 DPI
Anigen AIV Ag	1/8 (12.5)	3/8 (37.5)	8/8 (100)	8/8 (100)	8/8 (100)	8/8 (100)	0/8 (0.0)
Virus Isolation	4/8 (50.0)	5/8 (62.5)	8/8 (100)	8/8 (100)	8/8 (100)	8/8 (100)	4/8 (50.0)
RT-PCR	3/8 (37.5)	4/8 (50.0)	8/8 (100)	8/8 (100)	8/8 (100)	7/8 (87.5)	2/8 (25.0)

Anigen AIV Ag Test kit shows 100% correlation from 3-5 days post infection comparing Virus Isolation and RT-PCR

Evaluated at Rearch Unit of Green Cross Veterinary Products Co. Ltd., Korea

4. Comparison Study with other companies

Strain	EID50	Dil. virus	Anigen ICT*	"B" Flu Ag ICT	"I" AIV Ag ELISA
A/ty/lt/2962/03/H7N3 LPAI	108,5/ml	10 ⁻¹	+	+	+
		10 ⁻²	+	+	+
		10 ⁻³	+	+	+
		10 ⁻⁴	-	-	-
A/ty/lt/214845/02/H7N3 LPAI	107,5/ml	10 ⁻¹	+	+	+
		10 ⁻²	+	+	+
		10 ⁻³	+/-	+	+
		10 ⁻⁴	-	-	-
A/ty/lt/90302/05/H5N2 LPAI	107,5/ml	10 ⁻¹	+	+	+
		10 ⁻²	+	+	+
		10 ⁻³	+	+/-	+
		10 ⁻⁴	-	-	-
A/Chicken/Vietnam #8 (Y-irradiated) - A(H5N1)			+	NT**	NT
A/Duck/Vietnam #7A (Y-irradiated) - A(H5N1)			+	NT	NT
A/?/Vietnam #Pooled (Y-irradiated) - A(H5N1)			+	NT	NT

(*ICT: Immunochromatographic test, **: NT(Not Tested))

Matrix	Animal	Anigen ICT Kit	"B" Flu Ag ICT	"I" AIV Ag ELISA	Inoculation	Strain
Faeces	Duck	+	NT**	+/-	+	H11N9
Faeces	Duck	+	NT	+	+	H9N8
Tracheal	Turkey	+	+	+	+	H7N3
Tracheal	Turkey	+	+	+	+	H5N2
Tracheal	Turkey	+	+	+	+	H5N2
Faeces	Avian	+	NT	+	+	H5N3
Faeces	Turkey	-	NT	-	-	-
Tracheal	Avian	-	-	-	-	-
Lung	Pig	-	NT	-	-	-
Tracheal	Turkey	-	-	-	-	-

** NT (Not Tested)

"B" Flu Ag ICT: Not for animal use

5. Comparison of Rapid Diagnostic Test with AIV field samples

	Sensitivity (%)	Specificity (%)	Cohen's kappa
Anigen AIV Ag	76	97	0.71
A Company	70	90	0.58
B Company	63	97	0.57

*Evaluated by Australian Animal Health Laboratory (AAHL), Australian Quarantine and inspection Service (AQIS), Disease Investigation Centre Region IV, Yogyakarta, Indonesia

RAPID

Avian

H5 AIV Ag

Anigen Rapid H5 AIV Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of avian influenza type A subtype H5 virus antigen (High pathogenic avian influenza virus: HPAI) in avian cloaca swab, trachea swab, kidney swab and feces.



CAT.No.: RG15-05

Background

Avian influenza A viruses can be classified into low and high pathogenic forms based on the severity of the illness they cause. Avian influenza A subtype H5 and subtype H7 viruses can be distinguished as "high pathogenic" form on the basis of genetic features of the virus and the severity of the illness. Avian influenza A subtype H5 potentially has nine different subtypes (H5N1 ~ H5N9) and these subtypes can be highly pathogenic (HPAI) or low pathogenic (LPAI). The Influenza type A subtype H5 infection has been reported among humans and sometimes causing severe illness and death.


Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - Monoclonal anti-AIV H5 (Capture) - AIV H5 - Monoclonal anti-AIV H5 (Detector)
 - Monoclonal antibody against subtype H5 specific protein of avian influenza virus
- Purpose: Detection of Avian influenza type A subtype H5 virus (H5N1 ~ H5N9)
- Specimen: Avian cloca, trachea, kidney swab and feces
- Reading time: 5~10 minutes
- Sensitivity: 100% by farm (n=13), 76.6% by feces (n=115) vs. Virus isolation
- Specificity: 100% vs. Virus isolation, PCR (n=1,402)
- No cross reaction against NDV, IBV, ILTV, Pneumovirus, Reovirus, IBDV, MDV, and REV
- Detection limit: 104.9 EID50/ml or 0.125 HAU
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 25 Tests/Kit
- The Pigeons have very low possibility of AIV infection and the specimens from pigeon feces may show unusual result in the kit.


Special Features

- Optimal screening method of High pathogenic AIV subtype H5 in the farm or laboratory
- Detection of all avian influenza subtype H5 viruses (H5N1 ~ H5N9)
- High sensitivity compared with the virus isolation method
- High specificity: No cross-reaction with other avian viruses
- Detection of AIV H5 in Chickens, Ducks, Turkeys, and Geese
- World's first commercialized rapid test kit for detection of avian influenza virus subtype H5 since December 2003
- Easy test procedure: No technical expertise required, and Saving Labor & Time
- No additional equipment is required.
- Anigen AIV Ag test can be performed for initial screening, and Anigen H5 AIV Ag test for second line screening, HPAI subtype H5.

Test Procedures



In 5~10 minutes



1. Collect a cloaca feces swab, traceal swab, kidney swab or scattered wet feces swab.
2. Dilute the swab sample with assay diluent and drop 4 drops.
3. Interpret the test result.

Independent Comparison Study

Virus	Anigen H5 ICT	"BN" Flu A ICT
A/Chicken/Vietnam #8 (γ-irradiated) – A(H5N1)	Positive	Positive
A/Duck/Vietnam #7A (γ-irradiated) – A(H5N1)	Positive	Positive
A/?/Vietnam #Pooled (γ-irradiated) – A(H5N1)	Positive	Positive
A/Human/New Caledonia6/2002 – A(H3N2)	Negative	Positive
A/Human/Ulan Ude/1/2001- A(H3N2)	Negative	Positive

(*ICT: Immunochromatographic test)

Name of the farm infected with HPAI	Species	No. of random specimen	Anigen H5 ICT
			No. of Positive result
1.LKB	Egg Layers	3	2
2.KJT	Egg Layers	3	2
3.SKU	Egg Layers	3	1
4.SWS	Egg Layers	3	2
5.CMK	Egg Layers	5	5
6.STB	Egg Layers	6	6
7.KBK	Egg Layers	7	5
Total	Positive rate by Farms 7/7(100%) Positive rate by Specimens: 23/28(82.1%)		

Avian

AIV / H5 AIV Ag

Anigen Rapid AIV/H5 AIV Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of avian influenza type A viruses antigen and avian influenza subtype H5 virus antigen (High pathogenic avian influenza virus: HPAI) in avian cloaca swab, trachea swab, kidney swab or feces.



Background

Influenza viruses that infect birds are called "Avian Influenza virus (AIV)" with only influenza type A viruses infecting birds. Influenza type A viruses are divided into subtypes based on two proteins, hemagglutinin (HA) and neuraminidase (NA), on the surface of the virus. Avian influenza A viruses can be classified into low pathogenic and high pathogenic forms based on the severity of the illness they cause. The subtype H5 and subtype H7 viruses of Avian influenza A viruses can be distinguished as "high pathogenic" forms on the basis of genetic features of the virus and the severity of the illness. The AIV subtype H5 is potentially nine different subtypes (H5N1 ~ H5N9) and these subtypes can be highly pathogenic (HPAI) or low pathogenic (LPAI). The AIV subtype H5 infections have been reported among humans, sometimes causing severe illness and death.

Specification of AIV type A detection

- Principle: Immunochromatographic assay using Direct Sandwich Method
- Monoclonal anti-AIV (Capture) - AIV Type A - Monoclonal anti-AIV
- Monoclonal antibody against Type A specific nucleoprotein of Avian influenza virus
- Purpose: Detection of Avian influenza type A virus (H1~H16)
- Specimen: Avian cloaca, trachea, kidney or feces swab
- Sensitivity: 99.9% vs. HA, ELISA, 77.3% by feces vs. virus isolation
- Specificity: 99.9% vs. HA, PCR
- No cross reaction against NDV, IBV, ILTV, Pneumovirus, Reovirus, IBDV, MDV, REV
- Detection limit: 103.1 EID50/ml
- Reading time: 5~10 minutes


Specification of AIV subtype H5 detection

- Principle: Immunochromatographic assay using Direct Sandwich Method
- Monoclonal anti-AIV H5 (Capture) - AIV H5 - Monoclonal anti-AIV H5
- Monoclonal antibody against subtype H5 specific protein of avian influenza virus
- Purpose: Detection of Avian influenza type A subtype H5 virus (H5N1~H5N9)
- Specimen: Avian cloaca, trachea, kidney or feces swab
- Sensitivity: 100% by farm (n=13), 76.6% by feces vs. Virus isolation
- Specificity: 100% vs. Virus isolation, PCR
- No cross reaction against NDV, IBV, ILTV, Pneumovirus, Reovirus, IBDV, MDV, REV
- Detection limit: 104.5 EID50/ml or 0.06 HAU
- Reading time: 5~10 minutes

Special Features

- Optimal screening method of AIV Type A and High pathogenic subtype H5 in the field, laboratory
- High Sensitivity comparing with the virus isolation method
- High specificity: No cross-reaction with other avian viruses.
- Detection of AIV A and H5 in Chicken, Duck, Turkey, Goose
- World first commercialized rapid test kit for detection of avian influenza virus A and H5
- Easy test procedure: No technical expertise required and Labor & Time saving. No additional equipment is required.
- Fast Result: Result comes out within 5~10 minutes.
- Simultaneously AIV screening and HP AI subtype H5
- Shelf life: 24 months, -Storage temperature: 2~30°C
- Catalogue No.: RG 15-06

Test Procedures



Interpretation of the test result

1) Negative Result



2) AIV Type A Positive Result



3) AIV H5 Positive Result



4) Invalid Result



Anigen Rapid AIV Ab Test kit is a solid phase chromatographic immunoassay for the qualitative detection of antibody against Avian Influenza A in serum or egg yolk.



CAT.No.: RB25-01

Background

Wild birds and some poultry infected by HPAI show no clinical signs, and their virus excretion is too to be detected by an antigen capture immunoassay. To detect these latent AIV carriers, BIONOTE has developed the Anigen AIV Ab Rapid Test Kit as another for detecting AIV antibodies in these bird populations.

Specifications

- Principle: Competitive chromatographic immunoassay
- Purpose: Detection of AIV type A antibodies
- Species: Chicken, duck, turkey or wild birds
- Specimen: Serum or egg yolk.
- Reading time: 30 minutes
- Sensitivity: In case of HI 28, 96% (24/25)
- Specificity: In case of HI 20, 96.5% (55/57)
- No cross reaction against NDV, IBV, IBDV or MDV positive sera
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 10 Tests/kit

Special Features


- Optimal screening method of High pathogenic AIV in carrier birds
- No special equipment or training required
- Suitable test method at the lab. or farm.
- Applicable to various species
- The world's first antibody for AIV Rapid test Kit
- High sensitivity and specificity compared with the HI test

Sensitivity and Specificity

HI titer	Sensitivity		Specificity	
	ELISA Kit	Rapid Kit	ELISA Kit	Rapid Kit
Negative	2 ⁰		95.1% (77/81)	96.5% (55/57)
	2 ¹		100% (5/5)	100% (2/2)
	2 ²		38.5% (5/13)	87.5% (7/8)
Positive	2 ³	96.6% (28/29)	60.7% (17/28)	
	2 ⁴	98.6% (68/69)	88.0% (59/67)	
	2 ⁵	98% (96/98)	91.3% (84/92)	
	2 ⁶	100% (65/65)	95.2% (60/63)	
	2 ⁷	97.9% (46/47)	100% (43/43)	
	2 ⁸	100% (28/28)	96% (24/25)	

Test Procedures

- 1) Take a sample of poultry serum with a disposable dropper.
- 2) Slowly add one drop of serum to the sample hole and wait for 10 minutes, and then add 3 drops of assay diluent.
- 3) Test results will appear as purple band(s) in the result window of the test device.
- 4) Interpret test results at 20 minutes after assay diluent dropping.



* For egg yolk samples, 2ml of egg yolk should be mixed with 2ml of PBS (pH7.2) using a vortex agitator. And then the mixed egg yolk should be centrifuged at 3,000 rpm for 30 minutes. The supernatant is then used as the sample.

Avian NDV Ag

Anigen Rapid NDV Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of Newcastle Disease virus antigen in avian cloaca or trachea swab.



CAT.No.: RG15-03

Background

Newcastle Disease (ND) is characterized by sneezing, coughing, and nervous behavior. Affected birds may show tremors, circling, falling, twisting of the head and neck, or complete paralysis. Mortality reaches 90% in very young birds but adult mortality is very low. Among affected laying hens, egg quantity and quality drop sharply at first but usually return to former levels within four to eight weeks. ND can cause death even in vaccinated poultry and is characterised by its sudden appearance and rapid spread to nearby flocks. Diagnosis is usually made by virus isolation, serology and clinical signs. The hemagglutination inhibition (HI) test is used widely in ND virus serology, with cross-reactions in HI tests between ND virus and some of the other paramyxoviruses being reported. The newly developed Anigen Rapid NDV Ag test kit is highly specific to NDV.

Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
- Monoclonal anti-NDV (Capture) – NDV – Monoclonal anti-NDV (Detector)
- Purpose: All new NDV strains (Lasota, B1, VG/GA)
- Specimens: Avian cloaca or trachea swab
- Reading time: 5-10 minutes
- Sensitivity: 94.7% vs. RT-PCR
- Specificity: 96.4% vs. RT-PCR
- No cross reaction against AIV, IBV (infectious bronchitis), IBDV, Mycoplasma
- Detection limit: 0.125 HAU (Hemagglutination Unit)
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 10 Tests/Kit

Special Features

- NDV control is aided by fast and accurate diagnosis.
- Detection of all strains of NDV.
- World first commercialized rapid test kit for detection of NDV.
- Sample swabs are included for convenient test specimen sampling.
- Rapid test result: 10 minutes.
- Easy test procedure: No technical expertise required, Labor & Time saving.
- No additional equipment is required.

Test Procedures



Wait for
10 min.

NDV Ag Positive



NDV Ag Negative



- Collect a cloaca or trachea swab
- Dilute the swab sample with assay diluent and drop 4 drops
- Wait for 10 minutes and read the test result

Sensitivity & Specificity Study

Table 1) Sensitivity study compared with a commercial NDV RT-PCR Kit

	No. of Specimen	Anigen Rapid NDV Ag	RT-PCR
Positive result	18	+	+
	1	-	+
Sensitivity	19	18/19(94.7%)	19/19

Table 2) Specificity study compared with a commercial NDV RT-PCR Kit

	No. of Specimen	Anigen Rapid NDV	RT-PCR
Negative result	81	-	-
	3	+	-
Specificity	84	81/84(96.4%)	84/84

RAPID

Avian

IBDV Ag

Anigen Rapid IBDV Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of Infectious Bursal Disease (Gumboro) virus antigen in avian Bursa of Fabricius or feces.



Background

Infectious Bursal Disease (IBD), or Gumboro Disease, is a viral disease usually affecting young chickens 3 to 6 weeks old, and is spread through contaminated feed and water. The target organ of the virus is the Bursa of Fabricius, an important organ in the young chickens developing immune system. IBDV serotype 1 causes clinical disease in chickens younger than 10 weeks, with older chickens usually showing no clinical signs. IBDV serotype 2 is very widespread in turkeys and is sometimes found in chickens and ducks. In practice, a diagnosis can be indicated by the sudden onset of mortality in chickens between 2 and 8 weeks of age, and the presence of distinctive lesions in the Bursa of Fabricius and accompanying blood spots in the musculature of the breast and thigh of affected chickens

Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - Monoclonal anti-IBDV (Capture) – IBDV – Monoclonal anti-IBDV (Detector)
- Purpose: Detection of Infectious Bursal Disease virus antigen
- Specimen: Avian Bursa of Fabricius or cloaca swab
- Reading time: 5-10 minutes
- Sensitivity: 99.9% vs. RT-PCR
- Specificity: 96.6% vs. RT-PCR
- No cross reaction with AIV, IBV (Infectious bronchitis), NDV, Mycoplasma
- Detection limit: 104.0 EID50/ml
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 10 Tests/Kit

Special Features

- IBDV control is aided by fast and accurate diagnosis
- Good correlation with RT-PCR results
- Detection of all IBDV
- World first commercialized rapid test kit for detection of IBDV
- Convenient sampling of test specimens with a swab
- Rapid test result: 10 minutes.
- Easy test procedure: No technical expertise required, Saving time and labor
- No additional equipment is required

Sensitivity and Specificity Study

(1) Sensitivity study in 20 feces originated from IBDV positive chickens

Result	Anigen Rapid IBDV	RT-PCR
Positive	17	20
Negative	3	0
Sensitivity	17/20(85%)	20/20

(2) Sensitivity study in 12 Bursa of Fabricius originated from IBDV positive chickens


Result	Anigen Rapid IBDV	RT-PCR
Positive	12	12
Negative	0	0
Sensitivity	12/12(100%)	12/12

(3) Specificity study in negative feces from healthy chickens


Result	Anigen Rapid IBDV	RT-PCR
Positive	15	150
Negative	145	150
Specificity	145/150(96.6%)	150/150

Test Procedures

1) Swab method of Bursa of Fabricius



2) Cloacal Swab Method



3) Result Interpretation



RAPID

Bovine

B. brucella Ab

Anigen Rapid B. Brucella Ab Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of *Brucella abortus* antibody in whole blood, plasma, serum or milk.



CAT.No.: RB23-01

Background

Bovine brucellosis is commonly caused by *Brucella abortus* and less frequently by *B. melitensis*, and rarely by *B. suis*. Infection is widespread globally. Humans may be infected by contact with animals or animal products contaminated with these bacteria. Available serological tests include the Rose Bengal, ELISA, complement fixation test and tube agglutination test. However, these tests do not provide rapid diagnosis and each requires specialised laboratories and equipment. The new Anigen immunochromatographic rapid assay was developed to provide accurate, rapid and easy testing of *B. brucellosis*.


Specifications

- Principle: Immunochromatographic assay
 - (LPS capture)-(Antibody in blood, milk)-(Protein G detector)
 - LPS (Lipopolysaccharide) antigen originated from *B. abortus*
- Purpose: Detection of *Brucella abortus*, *melitensis* or *suis* antibodies
- Specimen: Bovine whole blood, serum, plasma and raw milk
- Reading time: 20 minutes
- Sensitivity: 97.9% vs. a commercial ELISA, 96.2% vs. Milk Ring Test.
- Specificity: 96.8% vs. a commercial ELISA, 89.7% vs. Milk Ring Test
- No cross reaction against *Brucella canis*, *Yersinia Enterocolitica*
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 30 multi-devices/Kit, 10 tests/Kit

Special Features

- Optimal screening method for Bovine brucellosis in the farm or laboratory
- Detection of antibodies against *Brucella abortus*, *melitensis*, *suis*
- High sensitivity versus a current confirmatory ELISA kit in blood, milk specimens
- Detection of *B. brucellosis* in goats, cows and oxen
- World's first commercialized rapid test kit for detection of *B. brucellosis*
- Easy test procedure: No technical expertise required
- Labor & Time saving with results in 20 minutes
- All components are included in the test kit. No additional equipment is required.
- Fast diagnosis helps to minimise losses to farmers and their cattle

Test Procedures



	Rose Bengal & Milk Ring Test	Anigen Rapid Test	Anigen's benefit
Specimen	-Rose Bengal: Serum -Milk Ring Test: Milk	Whole blood, Serum, Plasma, Raw milk available	All specimens are available
Testing Method	- Rose Bengal 1)Preparation of serum 2)1 drop of reagent of Rose Bengal + 1drop of serum : Mixing 3)Interpret the result within 4 min. - Milk Ring Test (MRT) 1)Adding of 1ml of milk into a tube 2)Adding of 30ul of MRT reagent 3)Waiting for 1 hour at 37 °C 4)Reading the result	1) Adding 1 drop of specimen 2) Adding one drop of assay diluent 3) Reading the test result at 20 min.	Easy test procedure
Reading Method	-Rose Bengal : Agglutination -MRT : Changing of blue color of Milk	Positive : 2 lines Negative : 1 line	Simple
Testing time	one day	20 Min.	Rapid Testing

Independent Comparison Study

1) Sensitivity and Specificity in serum

No. of Specimen		VLA B. Brucella ELISA, UK	
		Positive result	Negative result
Anigen Rapid B. Brucella Ab Test	Positive result	30	4
	Negative result	2	224
	Total	32	228
Sensitivity		30/32=93.8%	
Specificity		224/228=98.2%	

2) Sensitivity and Specificity in raw milk

Raw milk		"S" ELISA for B.Brucella Ab	
		Positive result	Negative result
Anigen Rapid B. Brucella Ab Test	Positive result	30	0
	Negative result	1	35
	Total	31	35
Sensitivity		30/31=96.7%	
Specificity		35/35=100%	

Raw milk		Milk Ring Test	
		Positive result	Negative result
Anigen Rapid B. Brucella Ab Test	Positive result	26	4
	Negative result	1	35
	Total	27	39
Sensitivity		26/27=96.2%	
Specificity		35/39=89.7%	

RAPID

Bovine

Bovine TB Ab

Anigen Rapid Bovine TB Ab Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of *Mycobacterium bovis* antibody in serum or plasma



CAT.No.: RB23-02

Background

Bovine tuberculosis is a chronic infectious disease caused by *Mycobacterium bovis*. Bovine tuberculosis infection in cattle is usually diagnosed in the live animal on the basis of delayed hypersensitivity reactions. Infection is often sub clinical. When present, clinical signs are not specifically distinctive of this disease. Generally the purified protein derivative (PPD) test has been known 65.6% sensitivity (Wood et al., Vet. Microbiol, 40:125-135. 1994) as non *M. bovis*, including *M. avium* and other acid-fast bacteria infect the cattle in nature. A more sensitive and accurate diagnostic tool has been required. Recently, the MPB 70 protein was revealed to be a highly species specific protein secreted by *Mycobacterium bovis*, Resulting in the development of the Anigen Rapid Bovine TB Ab test kit. The new diagnostic blood test using the MPB 70 protein is now available.

Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - (Rec. MPB70 Ag capture)-(Antibody in blood)-(Rec. MPB70 Ag detector)
 - Capture and detector material: Recombinant MPB70 antigen
- Purpose: Detection of antibodies (IgM, IgG) against *Mycobacterium bovis*
- Specimen: Bovine serum, plasma
- Reading time: 20 minutes
- Sensitivity: •90% vs. BTB positive confirmed by bacterial isolation
 - 85.1% vs. PPD test
- Specificity: 98.6% vs. PPD test
- No cross reaction against *Mycobacterium avium*, *Mycobacterium tuberculosis*
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 30 multi-devices/Kit, 10 individual devices/Kit

Special Features

- Optimal screening method of BTB in the farm or laboratory
- Higher sensitivity than the conventional PPD Test
- High specificity as it reacts only with the antibody against *M. bovis*
- No additional equipment is required
- Easy to use, saving time and labor

Clinical Evaluation Study

	Serum from naturally <i>M. bovis</i> infected cow (n=49)	Serum from artificially <i>M. bovis</i> inoculated cow n=20	Serum from <i>M. bovis</i> infected cow (n=69)
Positive rate of Anigen B.TB Ab	41/49 (83.7%)	18/20 (90%)	59/69 (85.5%)

What is the MPB70?

- MPB70 was identified as a major antigen from culture filtrate protein of *M. bovis*.
- MPB 70 is present predominantly in *M. bovis* and to a very limited degree in *M. tuberculosis* and other mycobacterial species (*M. avium* etc.).
- For that reason, Anigen BTB Ab test kit doesn't have cross-reaction with other species.


[Comparison Table of Anigen BTB Ab Rapid Test vs. PPD Test]

Virus	PPD Test	Anigen B. TB Ab Rapid Test	Anigen's Advantage
Principle	Delayed hypersensitivity with a purified protein derivatives	Chromatographic Immunoassay: Humoral immunity	
Specimen	Intradermal direct injection	Serum, plasma	
Method	1) Intradermal injection of bovine tuberculin PPD 2)The subsequent detection of swelling at the site of injection 3 days later	1) Adding 3 drop of specimen 2) Reading result at 20 min.	Easy Test procedure
Reading Method	By Observation of increases in skin-fold thickness - No more than 2 mm : Negative - 2 ~ 4mm : Inconclusive =>Another test after an interval of 42 days - More than 4mm : Positive	By Observation of line appearance 2 line => Positive 1 line => Negative	Simple
Testing time	Minimum 3 days	20 Min	Rapid Testing
Sensitivity	65% vs. Culture of <i>M. bovis</i> => High false positive rate due to the cross reaction with <i>M. avium</i> and others	90% vs. Culture of <i>M. bovis</i>	High Accurate
Early infection stage	False negative responses may occur	Less false negative responses may occur : By detection of IgM	Early detection
Late infection stage	False negative responses may occur	Less false negative responses may occur : By detection of IgG	

Test Procedures



Features of Test Result



RAPID

Bovine

FMD NSP Ab

Anigen Rapid FMD NSP Ab Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of Foot and Mouth Disease virus antibody in whole blood serum or plasma



CAT.No.: RB23-03

Background

Foot-and-Mouth (FMD) disease is a severe, highly contagious viral disease of cattle and swine. It also affects sheep, goats, deer, and other cloven-hooved ruminants. Signs of illness can appear after an incubation period of 1 to 8 days, but often develop within 3 days. There are FMD serotypes: O, A, C, SAT-1, SAT-2, SAT-3, and Asia-1. These subtypes show some regionality, and the O serotype is most common. The disease is characterized by fever and blister-like lesions followed by erosions on the tongue and lips in the mouth, on the teats, and between the hooves. It causes severe losses in the production of meat and milk.

Explanation of the test

The Anigen Rapid FMD NSP Ab Test Kit is an immunochromatographic assay for the qualitative detection of FMDV antibody in whole blood, plasma and serum from bovine, ovine and porcine animals. It can be used for identification of FMD field infection in cattle, pig, goat and sheep. The specially selected recombinant nonstructural protein (NSP) antigens are used as capture material in the test

Specifications

- Principle: Immunochromatographic assay
- Coated Antigen: Recombinant NSP antigen
- Specimens: Blood, Plasma, Serum
- Species: Cattle, Sheep, Goat, Pig
- Total Testing Time: 20~30 minutes
- Shelf life: 24 months
- Packing size: 10Tests/Kit
- Storage: 2~30°C

Advantages

- Mass screening in the farm or laboratory
- Confirmation of FMDV infection in non-vaccination herds
- Differential test of FMDV infected of vaccinated ones
- Excellent surveillance tool of FMDV
- Saving Labor & Time
- No additional equipment or technical expertise required

Sensitivity and Specificity study

1. Bovine

1) Specimen collected in outbreak region(Sensitivity)

		AniGen ELISA	Commercial ELISA
Anigen	+	217	270
Rapid	-	8	13
		225	283

2) Specimen collected in FMD free region(Specificity)

Anigen	+	507
Rapid	-	5
		512

3) Result

	Anigen Rapid FMD NSP	
	Sensitivity	Specificity
AniGen ELISA	94.4% (217/225)	99.0%
Commercial ELISA	95.4% (270/283)	

2. Swine

1) Specimen collected in outbreak region(Sensitivity)

		AniGen ELISA	Commercial ELISA
Anigen	+	20	20
Rapid	-	1	2
		21	22

2) Specimen collected in FMD free region(Specificity)


Anigen	+	566
Rapid	-	9
		575

3) Result

	Anigen Rapid FMD NSP	
	Sensitivity	Specificity
AniGen ELISA	95.2% (20/21)	98.4%
Commercial ELISA	90.9% (20/22)	


Test Procedures

(1) Serum & plasma specimen



- Add 5 ul of serum or plasma to the sample hole with a capillary tube
- Wait 1 minute
- Add 3 drops of the buffer
- Read the test result at 20~30 minutes

(2) Whole blood specimen



- Dispense 3 drops of the whole blood diluent into the dilution well
- Add 1 drop(25ul) of a whole blood sample with a disposable dropper
- Mix them 1 minute
- Extract the mixture, and add 5ul of the mixed sample to the sample hole with a capillary tube
- Wait 1 minute
- Dispense 3 drops of the buffer into the buffer well
- Read the test result at 20~30 minutes

RAPID

Bovine

Rota Ag

Anigen Rapid Rota Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of Rotavirus group A antigen in porcine, bovine, or canine feces



CAT.No.: RG18-03

Background

Rotavirus is the most common viral causative agent of diarrhea in pigs, cattle and dogs. Group A and B rotavirus are involved with group A being most prevalent and clinically important, containing several serotypes of differing virulence. All ages are susceptible. Dogs show somewhat opportunistic clinical signs. If neonatal pigs and cattle do not receive protective levels of maternal antibody, they are likely to develop profuse watery diarrhea in 12-48 hr. More commonly, the infection is endemic in a herd. Differential diagnosis of Rotavirus, TGEV and PEDV is important. Usually confirmatory diagnosis is based on histologic demonstration of villous atrophy in the jejunum or electron microscopy demonstration of virions in the intestinal contents. These methods are slow and require high priced specialized equipment and expertise. The Anigen Rapid Rota Ag test kit has been developed to provide fast, reliable, simple and sensitive detection of viral antigen in the intestinal mucosa or feces.

Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - Monoclonal anti-Rotavirus (Capture) – Rotavirus – Monoclonal anti-Rotavirus (Detector)
 - Monoclonal antibody against VP6 of Rotavirus
- Purpose: Detection of Rotavirus Group A antigen in pig, cattle and canine animal
- Specimen: Diarrheal feces
- Reading time: 10 minutes
- Sensitivity: 100% vs. RT-PCR
- Specificity: 98.5% vs. RT-PCR
- No cross reaction against PEDV, TGEV, PRC virus, BCV, E.coli, K88ab, K88ac, 987P, F41, Salmonella spp. Brachyspira hyodysenteriae, Clostridium perfringens, Coccidium, CPV, CCV
- Shelf life: 24 months
- Storage Temperature: 2~30°C
- Packing size: 10 individual devices/Kit

Independent Comparison Study

Clinical evaluation study of Anigen Rapid Rota Ag Test Kit

	Suckling pig	Calf
No. of specimen	12	18
No. of positive result of Antigen Rapid Rota Ag Test	5	6
No. of positive result of RT-PCR	5	6
Positive rate of Antigen Rapid Rota Ag Test	44%	33%
Positive result rate of RT-PCR	44%	33%
Agreement rate between Anigen Rota Ag Test and RT-PCR assay	100%	100%


Special Features

- Optimal detection method of Rotavirus group A antigen in the farm or laboratory.
- High sensitivity Vs RT-PCR
- Quick test result allows immediate preventive action to be taken.
- World's first commercialized rapid test kit for detection of Rotavirus antigen.
- Easy test procedure: No technical expert is required, saving labor & time
- All components are included in the test kit: Non-expert can perform testing.
- Fast diagnosis helps to decrease herd losses.

[Proper sample amount collected by a swab]



Test Procedures



RAPID

Swine

PED Ag

Anigen Rapid PED Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of Porcine Epidemic Diarrhea virus antigen in porcine feces.



CAT.No.: RG14-01

Background

Porcine epidemic diarrhea virus (PEDV) is a member of the coronavirus group which causes watery diarrhea, dehydration and high mortality in suckling pigs. Porcine Transmissible Gastroenteritis (TGE) virus is serologically unrelated with PEDV, but clinically these two virus infections are difficult to differentiate. The PEDV is transmitted by feces from infected pigs after oral uptake. Current available laboratory diagnosis in neonates is by direct immunofluorescence on cryostat sections of the small intestine or colon. ELISA and RT-PCR assays to detect viral antigens in feces or intestinal contents are useful, but require specialized instruments and laboratory personnel. The Anigen Rapid PED Ag Test Kit was developed for a non-expert to rapidly diagnose PEDV in farms or the laboratory.

Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - Monoclonal anti-PEDV (Capture) – PEDV – Monoclonal anti-PEDV (Detector)
 - Monoclonal antibody against Porcine Epidemic Diarrhea virus
- Purpose: Detection of Porcine Epidemic Diarrhea virus antigen in pig
- Specimen: Diarrhea feces
- Reading time: 10 minutes
- Sensitivity: 92% vs. RT-PCR.
- Specificity: 98% vs. RT-PCR
- No cross reaction against TGE virus, Rota virus, PRC virus, BCV, E.coli, K88ab, K88ac, 987P, F41, Salmonella spp. Brachyspira hyodysenteriae, Clostridium perfringens, Coccidium
- Shelf life: 24 months
- Detection limit: 1.0×10^3 TCID₅₀/ml or 0.5 HAU
- Storage Temperature: 2~30°C
- Packing size: 10 individual devices/Kit

Special Features

- It is the easiest way to confirm the cause of pig diarrhea
- Quick test results allow immediate preventive action to be taken
- High sensitivity gives confidence in results
- High specificity with no cross-reactivity
- No specialised equipment or personnel required
- One Step Testing - easy to use, Saving labor & time
- World's first & sole PED virus antigen Rapid Test Kit

Independent Clinical Evaluation Study


1. The Anigen Rapid PED Ag have tested with PED positive and negative clinical samples confirmed by RT-PCR at Korean National Veterinary Quarantine Service.

		RT-PCR		
		Positive	Negative	
Anigen Rapid PED Ag	Positive	164	4	Sensitivity: 92%
	Negative	15	184	Specificity: 98%

[Proper sample amount collected by a swab]



Test Procedures



RAPID

Swine

TGE Ag

Anigen Rapid TGE Ag Test Kit is a solid phase chromatographic immunoassay for the qualitative detection of Porcine Transmissible Gastroenteritis (TGE) Virus antigen in porcine diarrhea feces



CAT.No.: RG14-02

Background

Transmissible gastroenteritis (TGE) is a viral disease of the small intestine that causes vomiting and diarrhea in pigs of all ages. The infection spreads rapidly by aerosol or contact exposure. Severe epidemics are more common during winter due to survival of the virus in colder temperatures. Depending on the level of immunity and exposure, diarrhea may be mild in some litters but severe in others. Because TGE virus is easily spread during an epidemic by persons, animals, and other means, fast diagnosis and special care should be taken to prevent spread to unexposed groups of pigs and to neighboring herds. Current clinical signs in the epidemic form of TGE usually justify a presumptive diagnosis. In the mild endemic form, laboratory confirmation is required. The Anigen Rapid TGE Ag Test Kit provides fast and accurate confirmation of the presence of TGE viral antigens.


Specifications

- Principle: Immunochromatographic assay using Direct Sandwich Method
 - Monoclonal anti-TGEV (Capture) – TGEV – Monoclonal anti-TGEV(Detector)
 - Monoclonal antibody against S protein of Transmissible gastroenteritis virus
- Purpose: Detection of Transmissible gastroenteritis virus antigen in pig
- Specimen: Diarrhea feces
- Reading time: 10 minutes
- Sensitivity: 92.1% vs. RT-PCR & IFA
- Specificity: 95.2% vs. RT-PCR & IFA
- No cross reaction against PEDV, Rota virus, PRC virus, BCV, E.coli, K88ab, K88ac, 987P, F41, Salmonella spp. Brachyspira hyodysenteriae, Clostridium perfringens, Coccidium
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 10 individual devices/Kit

Special Features

- Optimal detection method of TGEV antigen in the farm or laboratory
- High sensitivity Vs current confirmatory RT-PCR
- Quick test result: Immediate preventive action is possible
- World's first commercialized rapid test kit for detection of TGEV antigen
- Easy test procedure: No technical expert is required, labor & time saving.
- All components are included in the test kit
- Fast diagnosis helps to decrease herd losses.

Test Procedures



Independent Comparison Study

1. Sensitivity Study: The Anigen Rapid TGE Ag Test Kit was tested with diarrhea feces to compare RT-PCR, IFA test result at Korean National Veterinary Quarantine Service

Pig Farm	No. of feces	Final diagnosis	No. of positive result	
			Anigen TGE Ag Kit	RT-PCR and IFA
A	12	TGEV	11	12
B	5	E.coli	0	0
C	12	TGEV	10	12
D	6	E.coli	0	0
E	14	TGEV	14	14
F	13	TGEV	11	13
G	15	TGEV + E.coli	14	15
H	9	TGEV	9	9
I	3	TGEV + E.coli	3	3
J	11	TGEV + E.coli	10	11
Total	100		82	89

- Sensitivity: 82/89=92.1% vs. RT-PCR and IFA

2. Specificity Study: The Anigen Rapid TGE Ag Test Kit was tested with diarrhea feces to compare RT-PCR, IFA test result at Korean National Veterinary Quarantine Service

Age of pig	No. of specimens	No. of negative results	
		Anigen TGE Ag	RT-PCR and IFA
Suckling pig	173	165	173
Growing pig	37	35	37
Total	210	200	210

- Specificity: 200/210=95.2% vs. RT-PCR and IFA

RAPID

Swine

TGE/PED Ag

Anigen Rapid TGE/PED Ag Test Kit is a solid phase chromatographic immunoassay for the simultaneously qualitative detection of Transmissible Gastroenteritis (TGE) virus antigen and Porcine Epidemic Diarrhea (PED) Virus antigen in porcine diarrhea feces.



Background

TGE and PED have a similar pathogenesis. The TGE and PED virus are excreted in feces. An acute outbreak of PED resembles the TGE outbreak. Clinical differentiation between TGE and PED is very difficult. However, as compared with TGE, the incubation period is longer (3-4 days), not all the litters of suckling pigs may become sick, and mortality in neonatal pigs is lower (average 50%).

TGE in its typical epidemic form causes a rapidly spreading diarrhea in animals of all ages with high mortality in neonates. With PED, the diarrhea spreads at a slower rate, and although diarrhea is seen in most of the litters, some litters may remain healthy even in the absence of immunity. Morbidity is 100% in older pigs. The Anigen Rapid TGE/PED Ag Test Kit detects viral antigens in feces or intestinal contents and easily differentiates between the PEDV and TGEV in pigs of all ages.

Specifications of TGE Ag

- Principle: Immunochromatographic assay using Direct Sandwich Method
- Monoclonal anti-TGEV (Capture) – TGEV – Monoclonal anti-TGEV (Detector)
- Monoclonal antibody against Transmissible Gastroenteritis virus
- Purpose: Detection of Porcine Transmissible Gastroenteritis virus antigen in pig
- Specimen: Diarrheal feces
- Reading time: 10 minutes
- Sensitivity: 92.1% vs. RT-PCR &IFA
- Specificity: 95.2% vs. RT-PCR &IFA
- No cross reaction against PEDV, Rota virus, PRC virus, BCV, E.coli, K88ab, K88ac, 987P, F41, Salmonella spp, Brachyspira hyodysenteriae, Clostridium perfringens, Coccidium
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 10 individual devices/Kit


Specifications of PED Ag

- Principle: Immunochromatographic assay using Direct Sandwich Method
- Monoclonal anti-PEDV(Capture) – PEDV – Monoclonal anti-PEDV(Detector)
- Monoclonal antibody against PED Virus
- Purpose: Detection of Porcine Epidemic Diarrhea virus antigen in pig
- Specimen: Diarrheal feces
- Reading time: 10 minutes
- Sensitivity: 92% vs. RT-PCR
- Specificity: 98% vs. RT-PCR
- No cross reaction against TGEV, Rota virus, PRC virus, BCV, E.coli, K88ab, K88ac, 987P, F41, Salmonella spp, Brachyspira hyodysenteriae, Clostridium perfringens, Coccidium
- Shelf life: 24 months
- Storage temperature: 2~30°C
- Packing size: 10 individual devices/Kit

Special Features

- World's first and sole rapid test kit for differential detection of TGEV and PEDV antigen.
- High sensitivity Vs current confirmatory RT-PCR
- Quick test result allows immediate preventive action to be taken.
- Easy test procedure: No technical expertise is required, saving labor & time
- Only one test sample is required for the differential testing.
- All components are included in the test kit
- Fast diagnosis helps reducing losses to herds.
- Anigen serves the first step to identify the virus causing the pig diarrhea

Test Procedures



Result Interpretation

1) TGE Ag and PED Ag Negative result	
2) Simultaneous TGE Ag and PED Ag Positive result	
3) TGE Ag Positive result	
4) PED Ag Positive result	

RAPID

AIV Ab ELISA

AniGen AIV Ab ELISA is a competitive Enzyme Linked Immunosorbent Assay for the qualitative detection of antibody against Avian Influenza A in serum, plasma and egg yolk.

Background

Wild birds and some poultry infected by HPAI show no clinical signs, and their virus excretion is too be detected by an antigen capture immunoassay. To detect these latent AIV carriers, BIONOTE has developed the AniGen AIV Ab ELISA as another for detecting AIV antibodies in these bird populations.

Specifications

- Principle: Competitive Enzyme Linked Immunosorbent Assay
[Nucleoprotein(NP) capture]-[AIV Ab in serum]-
[Monoclonal NP Ab-HRP detector]
- Purpose: Qualitative detection of AIV type A antibody
- Species: Chicken, duck, turkey, Geese, Swans, Quails, Guinea fowls
Grey partridges, Red partridges, Pheasants, Horses, pigs, Tern, Seal,
Gull, Mallard, Bar-headed Gull
- Specimen: Serum, plasma or egg yolk.
- Reading time: 45 minutes
- Sensitivity: Chicken 98.2%, Duck 97.5%, Turkey 97.1%
- Specificity: Chicken 97.3%, Duck 93.8%, Turkey 100%
- No cross reaction against NDV, IBV, IBDV or MDV positive sera
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features

- Optimal screening method of High pathogenic AIV in carrier birds
- Easy Test procedure : No sample dilution required
- Suitable test method at the lab.
- Applicable to various species
- Fast Test Result: Complete within 45 minutes
- High sensitivity and specificity compared with the HI test

Quick Procedure

1. Prepare AIV NP antigen coated test plate
2. Add 50ul of controls and sample to wells
3. Add 50ul of diluted enzyme conjugate to wells
4. Incubated plated for 30 minutes at 37°C
5. Wash plate 6 times
6. Add 100ul of substrate solution and incubated for 10 minutes at room temperature
7. Add 100ul of stop solution
8. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
9. PI value=[1-(OD sample/mean OD negative)] x 100

Sensitivity and Specificity

	Chicken	Duck	Turkey	Goose	Swan	Quail
Total No. of sera	1,613	178	213	25	4	46
PI value	50	50	85	50	50	51
Sensitivity (%)	98.2	97.5	97.1	100	100	100
Specificity (%)	97.3	93.8	100	100	100	97

	Guinea fowl	Grey Partridge	Red Partridge	Pheasant	Horse	Swine
Total No. of sera	19	38	5	18	63	266
PI value	19	50	50	50	50	50
Sensitivity (%)	87.5	100	Not tested	100	100	91.1
Specificity (%)	100	100	100	62.5	Not tested	97.4

*Validated of OIE AIV Reference Laboratory (IZS, Italy)



CAT.No.: EB45-02

Serotype Validation : This kit is valid in all subtype of AIV

H1N1	H2N2	H2N3	H3N2	H3N8	H4N6	H4N8	H5N1	H5N2	H5N3	H6N2
H7N1	H7N3	H7N7	H8N4	H9N2	H9N7	H10N1	H10N7	H11N6	H12N5	H13N6
H14N5	H15N6	H15N8	H16N3							

[Detection of AIV Ab in serum of Chicken that immunized purified inactivation virus]

Serum	Anigen Ab ELISA kit	HI*
PR/8/34 (H1N1)	+ (98)	102,400
Singapore/1/57 (H2N2)	+ (97)	16,384
Duck/Ukraine/1/63 (H3N8)	+ (99)	40,960
Duck/Czechoslovakia/56 (H4N6)	+ (99)	16,384
Tern/South Africa/61 (H5N3)	+ (99)	10,240
Turkey/Massachusetts/3740/65 (H6N2)	+ (94)	10,240
Seal/Massachusetts/1/80 (H7N7)	+ (98)	12,800
Turkey/Ontario/6118/68 (H8N4)	+ (99)	8,000
Turkey/Wisconsin/66 (H9N2)	+ (99)	16,000
Chicken/Germany/N/49 (H10N7)	+ (99)	51,200
Duck/England/1/56 (H11N6)	+ (99)	10,240
Duck/Alberta/60/76 (H12N5)	+ (99)	16,834
Gull/Maryland/704/77 (H13N6)	+ (84)	51,200
Mallard/Astrakhan/263/82 (H14N5)	+ (99)	10,240
Duck/Australia/341/83 (H15N8)	+ (98)	5,120
Bar-headed Gull/Sweden/5/99 (H16N3)	+ (99)	16,384
NDV Miyadera	- (-20)	131,072

*HI used each virus antigen that immunized

[Detection of AIV Ab in serum post experimental infection]

Serum	No	Anigen Ab ELISA kit				HI*		
		pre	1w	2w	3w	pre	1w	2w
Turkey/England/63 (H7N3) (HP AI)	1	- (-27)	+ (96)	+ (95)	+ (98)	<20	320	640
	2	- (-23)	+ (97)	+ (98)	+ (99)	<20	320	640
	6	- (-25)	+ (95)	+ (97)	+ (98)	<20	320	640
	7	- (-26)	+ (98)	+ (98)	+ (99)	<20	320	1280
	8	- (-22)	+ (90)	+ (96)	+ (98)	<20	320	640
	9	No test	+ (97)	+ (97)	+ (99)	<20	160	320
Chicken/Tbaraki/1/05 (H5N2) (LPAI)	11	No test	- (41)	+ (58)	+ (98)	<20	160	640
	12	No test	+ (58)	+ (69)	+ (99)	<20	160	640

*HI used each virus antigen that immunized

*Pre : Pre immunization, W : week(s)

*Validated of OIE AIV Reference Laboratory (Hokkaido Univ. Japan)

ELISA

Avian

H5 AIV Ab ELISA

AniGen H5 AIV Ab ELISA is a competitive Enzyme Linked Immunosorbent Assay for the qualitative detection of antibody against Avian Influenza A, H5 subtype in serum, plasma and egg yolk.



CAT.No.: EB 45-03

Background

Wild birds and some poultry infected by HPAI subtype H5 show no clinical signs, and their virus excretion is too be detected by an antigen capture immunoassay. To detect these latent H5 AIV carriers, BIONOTE has developed the AniGen H5 AIV Ab ELISA as another for detecting AIV subtype H5 antibodies in these bird populations.

Specifications

- Principle: Competitive Enzyme Linked Immunosorbent Assay [Hemagglutinin 5 capture]-[H5 AIV Ab in serum]- [Monoclonal Hemagglutinin 5 Ab-HRP detector]
- Purpose: Qualitative detection of AIV subtype H5 antibodies
- Species: Chicken, duck, Quail
- Specimen: Serum, plasma or egg yolk.
- Reading time: 45 minutes
- Sensitivity: Chicken 96.9% (Vs HI)
- Specificity: Chicken 100 % (Vs HI)
- No cross reaction against other AIV subtype antibodies
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features

- Optimal screening method of High pathogenic AIV, H5 in carrier birds and vaccination
- Easy Test procedure : No sample dilution required
- Suitable test method at the lab.
- Applicable to various species
- Fast Test Result: Complete within 45 minutes
- High sensitivity and specificity compared with the HI test Serotype Validation

Serotype Validation

Antibody	PI	Interpretation	Antibody	PI	Interpretation
H1N1	5	-	H5N2	95	+
H2N9	7	-	H5N3	98	+
H3N2	10	-	H7N7	11	-
H5N1	92	+	H9N2	13	-

Quick Procedure

1. Prepare AIV H5 antigen coated test plate
2. Add 50ul of controls and sample to wells
3. Add 50ul of diluted enzyme conjugate to wells
4. Incubated plated for 30 minutes at 37°C
5. Wash plate 6 times
6. Add 100ul of substrate solution and incubated for 10 minutes at room temperature
7. Add 100ul of stop solution
8. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
9. PI value=[1-(OD sample/mean OD negative)] x 100

ELISA

Avian

NDV velo Ab ELISA

AniGen NDV Velo Ab ELISA is an Indirect Enzyme Linked Immunosorbent Assay for the qualitative detection of antibodies against Velogenic Newcastle disease virus in avian serum or plasma



Background

Newcastle disease is a contagious bird disease affecting many domestic and wild avian species. NDV strains can be categorised as velogenic (highly virulent), mesogenic (intermediate virulence) or lentogenic (nonvirulent). Velogenic strains produce severe nervous and respiratory signs, spread rapidly and cause up to 90% mortality. Mesogenic strains cause coughing, affect egg quality and production and result in up to 10% mortality. Lentogenic strains produce mild signs with negligible mortality.

In acute cases, the death is very sudden, and, in the beginning of the outbreak, the remaining birds do not seem to be sick. In flock with good immunity, however, the signs (respiratory and digestive) are mild and progressive, and are followed after 7 days by nervous symptoms, especially twisted heads.

Specifications

- Principle: Indirect Enzyme Linked Immunosorbent Assay [Recombinant NDV antigen (Capture)]-[velogenic NDV antibodies in sample]-[Anti Avian IgG-HRP detector]
- Purpose: Qualitative detection of antibodies against Velogenic Newcastle disease virus
- Specimen: Serum or plasma
- Reading time: 105 minutes
- Sensitivity: 95.0%
- Specificity: 98 %
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features

- Differential test of Velogenic strain and other strains
- The first ELISA test kit for detection of antibodies against velogenic NDV virus
- Simple test procedure
- High sensitivity and specificity
- Fast Test result: Within 2 hours

Performance study

AniGen NDV Velo Ab ELISA	NDV Velogenic strain inoculation to vaccinated chickens					
	DPI 0	DPI 5	DPI 7	DPI 10	DPI 14	DPI 18
	0%	33%	70%	58%	25%	25%

* Test group: 12 chickens * Control group: 3 chickens

It is almost not possible to study sensitivity because there are no commercial tests or method confirmed for detection of the antibody by NDV strains.

Quick Procedure

1. Prepare NDV velo antigen coated test plate
2. Dilute test sample and control with sample diluents
3. Add 100ul of diluted sample and control to wells
4. Incubated plated for 60 minutes at room temperature.
5. Wash plate 5 times
6. Add 100ul of enzyme conjugate to wells
7. Add 100ul of substrate solution and incubated for 30 minutes at room temperature
8. Add 100ul of stop solution
9. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
10. $S/P \text{ value} = \frac{[OD \text{ sample} - \text{mean OD negative}]}{[\text{mean OD positive} - \text{mean OD negative}] \times 100}$

ELISA

Bovine

B.brucella Ab ELISA

AniGen B. Brucella Ab ELISA is an Indirect Enzyme Linked Immunosorbent Assay for the qualitative detection of Brucella abortus and melitenis antibody in bovine plasma, serum or milk.



CAT.No.: EB43-01

Background

Bovine brucellosis is a highly contagious bacterial disease, almost exclusively caused by Brucella abortus causing late term-abortion and infertility in cattle. This disease that affects many animal species is a zoonose of economic importance, as well as economic importance, as well as a public health hazard.

Specifications

- Principle: Indirect Enzyme Linked Immunosorbent Assay [LPS of Brucella abortus (Capture)]-[Brucella antibody in sample]-[Monoclonal Anti Bovine IgG-HRP detector]
- Purpose: Qualitative detection of Brucella abortus and melitenis
- Species: Bovine
- Specimen: plasma, serum or milk
- Reading time: 120 minutes
- Sensitivity: serum 100%, milk 100% (Vs Commercial ELISA)
- Specificity: Serum 97.9 %, milk 99.1% (Vs Commercial ELISA),
- No cross reaction against yersinia enterocolitica
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features

- High Sensitivity and Specificity in serum, plasma, or milk.
- Fully meets the requirement of EU Council Directive 64/432/EEC Annex C as last amended in March, 2002 and the OIE Manual of Diagnostics Tests and Vaccine Fifth Edition 2004 Chapter 2.3.1
- Standardized by OIE standard sera (OIESS and relevant OIE wp and OIEsp)
- Adapted to variety of automated robotic samplers and ELISA processor systems
- Individual or pool sample available

Quick Procedure

1. Prepare LPS of B. abortus coated test plate
2. Dilute the test samples 1:50 with sample diluents (Do not dilute controls and Milk)
3. Add 100ul of positive control, negative control, and diluted test samples.
4. Incubated plate for 60 minutes at 37°C
5. Wash plate 5 times
6. Add 100ul of diluted enzyme conjugate to wells
7. Incubated plate for 30minutes at 37°C
8. Wash plate 5 times
9. Add 100ul of substrate solution and incubated for 15 minutes at room temperature
10. Add 100ul of stop solution
11. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
12. %P value=[OD sample/mean OD positive)] x 100

Comparative Sensitivity and Specificity study

(1) Study in serum

Specimens : Total 271 of positive and negative cow serum from KNVRQS*

		AniGen ELISA			
			+	-	
Commercial ELISA	+	127	0	127	
	-	3	141	144	
		130	141		

Result

Sensitivity	100%
Specificity	97.9%

(2) Study in raw milk

Specimens : Total 271 of positive and negative cow milk from KNVRQS*

		AniGen ELISA			
			+	-	
Commercial ELISA	+	68	0	68	
	-	1	116	117	
		69	116		

Result

Sensitivity	100%
Specificity	99.1%

Bovine

BTB Ab ELISA

AniGen BTB Ab ELISA is a direct Enzyme Linked Immunosorbent Assay for the qualitative detection of *M. bovis* antibodies in bovine serum



CAT.No.: EB43-03

Background

Bovine tuberculosis infection in cattle is usually diagnosed in live animals on the basis of delayed hypersensitivity reactions (Intradermal skin test). The infection is often subclinical; when present, clinical signs are not specifically distinctive of this disease and might include weakness, anorexia, emaciation, dyspnoea, enlargement of lymph nodes, and cough, particularly with advanced tuberculosis. After death, it is diagnosed by post-mortem examination and histopathological and bacteriological techniques. DNA probe and polymerase chain reaction (PCR) techniques may also be used. Traditional bacterial culture diagnosis remains the routine method for confirmation of infection. The intradermal skin test is the method of choice for routine screening although it is time and labor intensive. The AniGen BTB Ab ELISA is an easy and fast serological diagnostic test intended for bovine tuberculosis mass screening before the intradermal skin test; the kit detects *M. bovis* antibody in bovine serum.

Specifications

- Principle: Direct Enzyme Linked Immunosorbent Assay [Recombinant *M. bovis* antigen (Capture)]-[*M. bovis* antibody in sample]-[Recombinant *M. bovis* antigen -HRP detector]
- Purpose: Qualitative detection of *M. bovis* antibody
- Species: Bovine
- Specimen: serum
- Reading time: 75 minutes
- Sensitivity: 87.3% (Vs isolation or intradermal skin test)
- Specificity: 98.9% (Vs isolation or intradermal skin test)
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features

- The world's first reliable BTB Ab ELISA for serological diagnosis
- Useful mass screening method for Bovine tuberculosis before intradermal skin test: Save time and labor cost
- Unique recombinant specific antigens provide high sensitivity and specificity
- Fast test result without sample dilution and secondary washing

Quick Procedure

1. Prepare BTB antigen coated test plate
2. Add 50ul of controls and sample to wells
3. Add 50ul of diluted enzyme conjugate to wells
4. Incubated plated for 60 minutes at 37°C
5. Wash plate 6 times
6. Add 100ul of substrate solution and incubated for 15 minutes at room temperature
7. Add 100ul of stop solution
8. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
9. $S/P \text{ value} = \frac{[OD \text{ sample} - \text{mean OD negative}]}{[\text{mean OD positive} - \text{mean negative}] \times 100}$

Diagnosis Flow

How can the AniGen BTB Ab ELISA be used effectively?

The individual Intradermal skin test is the standard screening method of diagnosis in the presence of live cattle. However, Intradermal skin test is more time and labor intensive compared to ELISA test. Therefore, it is recommended to use AniGen BTB Ab ELISA as a mass screening diagnostic method. ELISA positive cases are recommended to be tested again by Intradermal skin test.

1

AniGen BTB Ab ELISA test for mass screening

2

Intradermal skin test for suspected cattle

3

IFN- γ test for confirmation

ELISA



Bovine

BTB Ab ELISA

Performance data

Sensitivity

	AniGen BTB Ab ELISA			
	Positive	Negative	Total	Sensitivity
Positive serum confirmed by virus isolation	18	2	20	90.0%(18/20)
Positive serum confirmed by Intradermal skin test	113	17	130	86.9%(113/130)
Total	131	19	150	87.3%(131/150)

Specificity

	AniGen BTB Ab ELISA			
	Positive	Negative	Total	Specificity
Negative serum confirmed by epidemiological study	12	754	766	98.4%(754/766)
Negative serum confirmed by Intradermal skin test	80	8,220	8,300	99.0%(8,220/8,300)
Total	92	8,974	9,066	98.9%(8,974/9,066)

Q & A

1. When can the M. Bovis be detected by Intradermal skin test after M. Bovis infection?

Generally, the bovine tuberculosis can be detected by the intradermal skin test from 3~6 weeks after infection. The Intradermal skin test can detect M. Bovis at the early, middle, and late stage of a disease but not during the early stage.

2. When can the AniGen BTB Ab ELISA detect the anti- M. bovis antibody after M. bovis infection?

Generally, the bovine tuberculosis can be detected by the AniGen BTB Ab ELISA from 4~8.5 weeks after infection. The AniGen BTB ELISA can detect at the middle, late stage, and the last stage of a disease but not during the early stage.

3. What are the reasons for possible discrepancies between AniGen BTB Ab ELISA result and Intradermal skin test result?

- 1) In the case of a positive intradermal skin test, and a negative ELISA results.
 - If the cattle is truly infected, the early stage of M. bovis infection is stronger in cellular immunity, while it is weak in humoral immunity. Therefore, it is possible that Intradermal skin test that detects cellular immunity is positive and ELISA which detects humoral immunity is negative.
 - If the cattle is not infected, the false positive result is due to a non-specific response of the Intradermal skin test.
- 2) In the case of a negative Intradermal skin test, and a positive ELISA results.
 - If the cattle is truly infected, the last stage of M. bovis infection is stronger in humoral immunity, while it is weak in cellular immunity. Therefore, it is possible that Intradermal skin test that detects cellular immunity is negative and ELISA that detects humoral immunity is positive.
 - If the cattle is not infected, the false positive result is due to a non-specific response of ELISA.

Bovine

FMD NSP Ab ELISA

AniGen FMD NSP Ab ELISA is a Competitive Enzyme Linked Immunosorbent Assay for the qualitative detection of FMD NSP antibodies in Cattle, Sheep, Goat, Pig serum or plasma.



Background

Foot-and-mouth disease (FMD) is a highly contagious viral infection primarily of cloven-hoofed domestic animals, such as cattle, pigs, sheep, goats, deer, and water buffalo. The FMD virus has seven serotypes : A, O, C, Asia 1, and Southern African Territories (SAT) 1,2 and 3.

In many countries the disease is controlled by vaccinations that consist of (partly) purified structural proteins (SP) of the FMD virus, and therefore vaccinated animals only elicit antibodies directed against the structural proteins.

Non structural protein (NSP) is expressed only by replicating viruses, and inactivated vaccines are purified to remove the cellular proteins and NSP. Therefore only animals that have been infected with live virus should develop antibodies against NSP.

The detection of antibody to the NSP is the single most reliable indicator of FMDV infection (Lubroth, 1995, De Diego, 1997)

It is known that vaccinated animals, which are exposed to infection, can become persistently infected with FMDV without ever showing clinical signs (Mackay, D.K.J., et. Al., 1998)

In countries that use vaccination to control FMDV outbreaks, it is important to differentiate between antibodies against SP and NSP in to discriminate between infections from the field and immune response to vaccination.

Specifications

- Principle: Competitive Enzyme Linked Immunosorbent Assay [Recombinant FMD 3 ABC antigen (Capture)]-[FMD NSP antigen in sample]-[Monoclonal anti-FMD NSP antibody -HRP detector]
- Purpose: Qualitative detection of FMD NSP antibodies
- Species: Cattle, Sheep, Goat, Pig
- Specimen: serum or plasma
- Reading time: 120 minutes
- Sensitivity: Cattle 93.8%
- Specificity: Cattle 99.9%, Pig 99.9%, Sheep & goat 100%
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features



- Differential test of FMD infected or vaccinated
- High accuracy equivalent to a world standard ELISA kit
- Easy test procedure: No serum pre-dilution required
- Cost effective: No requirement for an uncoated microplate for serum pre-dilution
- Fast Test result: Within 2 hours

Quick Procedure

1. Prepare FMD NSP antigen coated test plate
2. Add 50ul of controls and sample to wells
3. Add 50ul of diluted enzyme conjugate to wells
4. Incubated plated for 90 minutes at 37°C
5. Wash plate 6 times
6. Add 100ul of substrate solution and incubated for 15 minutes at room temperature
7. Add 100ul of stop solution
8. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
9. PI value=[1-(OD sample/mean OD negative)] x 100

ELISA

FMD Immune response of infected & vaccinated cattle



Bovine

FMD NSP Ab ELISA

Differentiation Infected and Vaccinated.

FMDV Status	Serological Test result	
	LPB,VN Test	NSP test
Vaccinated, infected	Positive	Positive
Vaccinated, Not infected	Positive	Negative
Non vaccinated, Infected	Positive	Positive
Non vaccinated, Not infected	Negative	Negative

Performance characteristics

1) Cross reactivity to each FMDV serotype

Serotype	Detection Rate (positive no./ tested no.)
Asia 1	2/2
SAT 1	3/3
SAT 2	5/5
SAT 3	2/2
Type A	13/13
Type C	3/3
Type O	6/6

AniGen FMD NSP Ab ELISA detects antibodies against all FMDV serotypes
(The performance evaluation was performed in OIE FMD Ref. Laboratory)

2) Comparative Sensitivity study in experimentally infected animals

Species	Post Infection	Detection Rate (positive no./ tested no.)		
		AniGen	Company I	Company A
Cattle	5days	0/17	-	0/17
	7days	70/92(76.1%)	18/47 (38.3%)	54/92 (58.7%)
	14~26 days	30/32(93.8%)	30/32(93.8%)	31/32 (96.9%)
Sheep and Goat	7 months	3/6	-	3/6
	10-12 days	5/5	2/5	5/5
		0/3*	0/3*	0/3*
	21 days	3/3	2/3	3/3
		1/1*	1/1*	1/1*
Pig	7months	7/10	-	7/10
	10-12 days	4/5	-	-

AniGen FMD NSP Ab ELISA detects antibodies against nonstructural protein from 7 days to 7 months after infection
(*Experimentally contact challenge animal group, The performance evaluation was performed in OIE FMD Ref. Laboratory)

ELISA

Swine

CSFV Ab ELISA

AniGen CSFV Ab ELISA is an Indirect Enzyme Linked Immunosorbent Assay for the qualitative detection of antibodies against classical swine fever virus in swine serum



CAT.No.: EB44-13

Background

The infectious agent responsible is a virus CSFV (previously called hog cholera virus) of the genus Pestivirus in the family Flaviviridae. CSFV is closely related to the ruminant pestiviruses which cause Bovine Viral Diarrhoea (BVDV) and Border Disease (BDV). The effect of different CSFV strains varies widely, leading to a wide range of clinical signs. Highly virulent strains correlate with acute, obvious disease and high mortality, including neurological signs and hemorrhages within the skin. Less virulent strains can give rise to subacute or chronic infections that may escape detection, while still causing abortions and stillbirths. In these cases herds in high risk areas are usually serologically tested on a thorough statistical basis. Infected piglets born to infected but subclinical sows help maintain the disease within a population. Other signs can include lethargy, fever, immunosuppression, chronic diarrhea and secondary respiratory infections.

Specifications

- Principle: Indirect Enzyme Linked Immunosorbent Assay [Recombinant CSFV antigen (Capture)]-[CSFV antibodies in sample]-[Anti swine IgG-HRP detector]
- Purpose: Qualitative detection of antibodies against classical swine fever virus
- Species: Swine
- Specimen: Serum
- Reading time: 75 minutes
- Sensitivity: 99.3% (Vs Commercial ELISA)
- Specificity: 99.7 % (Vs Commercial ELISA)
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Tests/kit, 960 Test/ kit

Special Features

- No cross reaction with other swine disease positive sera.
- Simple test procedure
- High sensitivity and specificity
- Fast Test result: Within 2 hours

Quick Procedure

1. Prepare CSFV antigen coated test plate
2. Dilute test sample and control with sample diluents
3. Add 100ul of diluted sample and control to wells
4. Incubated plated for 30 minutes at room temperature.
5. Wash plate 5 times
6. Add 100ul of enzyme conjugate to wells
7. Wash plate 5 times
8. Add 100ul of substrate solution and incubated for 15 minutes at room temperature
9. Add 100ul of stop solution
10. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
11. $S/P \text{ value} = \frac{[OD \text{ sample} - \text{mean OD negative}]}{[\text{mean OD positive} - \text{mean OD negative}] \times 100}$

ELISA

Sensitivity and specificity

		Commercial ELISA		
		Positive	Negative	Total
AniGen CSFV Ab ELISA	Positive	156	0	156
	Negative	1	210	211
	Total	157	210	367

AniGen CSFV Ab ELISA has high Sensitivity 99.3% and high Specificity 99.7

Swine

PRRS Ab ELISA 4.0

AniGen PRRS Ab ELISA 4.0 is an Indirect Enzyme Linked Immunosorbent Assay for the qualitative detection of antibodies against PRRSV(EU, NA, and Korean strain) in swine serum.

Background

Porcine reproductive and respiratory syndrome (PRRS) is characterized by reproductive failure of sows and respiratory problems of piglets and growing pigs. The disease is caused by the PRRS virus, a virus currently classified as a member of the order Nidovirales, family Arteriviridae, genus Arterivirus. The primary target cell of the virus is the alveolar macrophage of the pig. Two major antigenic types of the virus exist, the European(EU) and the North American(NA) strain. The virus is primarily transmitted via infected pigs but also by feces, urine, semen and fomites.

Specifications

- Principle: Indirect Enzyme Linked Immunosorbent Assay [Recombinant PRRSV antigen (Capture)]-[PRRS Ab in sample]-[Polyclonal anti-swine IgG -HRP detector]
- Purpose: Qualitative detection of antibodies against PRRSV(EU, NA, and Korean strain) -Species: Cattle, Sheep, Goat, Pig
- Specimen: serum
- Reading time: 75 minutes
- Sensitivity: 98.7 (Vs IFA)
- Specificity: 99.7 (Vs IFA)
- Shelf life: 12 months
- Storage temperature: 2~8°C
- Packing size: 96 Tests/kit, 480 Test/kits, 960 Test/ kit

Special Features

- Antibodies against European strain, North America strain and Korean Strain can be detected.
- No cross reaction with other swine disease positive sera.
- Simple test procedure
- High sensitivity and specificity
- No more singleton reactor (false positive)
- Fast Test result: Within 2 hours

Quick Procedure

1. Prepare PRRSV antigen coated test plate
2. Dilute test sample and controls with sample diluents
3. Add 100ul of diluted sample and controls to wells
4. Incubated plated for 30 minutes at room temperature.
5. Wash plate 5 times
6. Add 100ul of enzyme conjugate to wells
7. Wash plate 5 times
8. Add 100ul of substrate solution and incubated for 15 minutes at room temperature
9. Add 100ul of stop solution
10. Measure the optical density (OD) at 450 nm with reference wavelength at 620nm
11. $S/P \text{ value} = \frac{[OD \text{ sample} - \text{mean OD negative}]}{[\text{mean OD positive} - \text{mean OD negative}] \times 100}$

PRRS strain detectability study

1) PRRS European strain

	Farm	No.of Samples	AniGen	Commercial Kit A
			%	%
Sensitivity	EU-1	10	100% (10/10)	90% (9/10)
	EU-2	6	100% (6/6)	100% (6/6)
	EU-3	9	100% (9/9)	100% (9/9)
	EU-4	8	87.5% (7/8)	100% (8/8)
	EU-5	8	100% (8/8)	75% (6/8)
	EU-6	28	64.3% (18/28)	75% (21/28)
	EU-7	29	100% (29/29)	96% (28/29)
Specificity	Farm A ¹⁾	14	93% (13/14)	93% (13/14)
	Multiple farms ²⁾	17	100% (17/17)	23.5% (4/17)

1) Farm A : 14 Samples from 1 farm show singleton reactor.

2) Multiple Farms : 17 Samples from several farms show singleton reactor.

2) PRRS Korean strain

Inoculation strain	IFA	DPI 3 weeks		
		AniGen PRRS Ab ELISA 4.0		Results
		S/P		
VR2332 like virus	No.1	1:64	1.5	Positive
	No.2	1:64	1.4	Positive
NSP2-deleted mutant		1:64	1.0	Positive
NA mutant (ORF5 mutant, Prevalent in Korea)	No.1	1:64	1.2	Positive
	No.2	1:64	1.3	Positive
EU Lelystad like virus		1:64	0.8	Positive

ELISA



CAT.No.: EB44-04

Study on PRRS standard antisera

Standard Sera	IFA		AniGen PRRS Ab ELISA 4.0			Commercial Kit A		
	Titer	Result	OD	S/P	Result	OD	S/P	Result
Standard antiserum (European strain)	>256	Positive	3.0060	4.8	Positive	0.8830	2.1	Positive
Standard antiserum (North American strain)	>256	Positive	2.2350	3.5	Positive	1.0270	2.4	Positive
Negative control serum	<20	Negative	0.1550	0.0	Negative	0.1440	0.0	Negative

AniGen PRRS Ab ELISA 4.0 can detect standard antiserum from both European and North American strain. *S/P ratio \geq 0.4: positive

Sensitivity and specificity study



		IFA		
		Positive	Negative	Total
AniGen PRRS Ab ELISA 4.0	Positive	324	1	325
	Negative	4	333	337
	Total	328	334	662

AniGen PRRS Ab ELISA 4.0 has high Sensitivity, 98.7% and high Specificity, 99.7%

Seroconversion study after challenge of PRRSV

1) EU strain Antibody response after challenge with PRRS European strain

(Sample No. EU-01, EU-02)




* DPI: Day Post Inoculation

AniGen PRRS Ab ELISA 4.0 can detect PRRS Ab from DPI 7 with PRRS European strain.

2) NA strain Antibody response after challenge with PRRS North American strain

(Sample No. NA-01, NA-02)



AniGen PRRS Ab ELISA 4.0 can detect PRRS Ab from DPI 7 with North American strain.

Avian General PCR Kits AIV detection Kit



PCR

Background

Avian Influenza virus are divided into subtypes based on two proteins, hemagglutinin(HA) and neuraminidase(NA), on the surface of the virus. AIV is spread primarily through direct form infected birds to healthy birds, and through indirect contact with contact with contaminated equipment and materials. The virus is excreted through infected birds' feces and secretions from their noses, mouths, and eyes. BIONOTE AIV Detection Test Kit is a complete, easy PCR kit for detection of Sepecific protein gene from Avian influenza type A virus (H1~H16 genes) BIONOTE AIV detection Kit is very accurate and high efficient, and RT-reaction can be done from 50fg to 500ng of template RNA. It is developed with the best condition of synthesis first-strand cDNA, so it is useful for check a low copy of DNA transcription.


Specifications

- ✓ **High Sensitivity and Specificity** : BIONOTE owning optimal& unique primers and probes system guarantees a highly specific result.
- ✓ **Excellent reproducibility**
- ✓ **Easy to use** : BIONOTE AIV Detection Kits contain all components for the detection of AIV.
- ✓ **Saving labor costs& time in the laboratory**
- ✓ **Includes positive and negative control**

Materials in AIV Detection kit (96tests/Kits)

Contents	Volume
Negative Control	0.5ml
Positive Control	0.5ml
20X RT Mix	0.1ml
20X PCR Enzyme Mix	0.1ml
Detection Solution	0.8ml

Interpretation of the AIV detection Kit



(S1: Negative sample, S2: Positive sample, BL: blank, P/C: Positive control, N/C: Negative control), 100b: 100bp ladder)

Type of AIV General PCR Kit

Product	Diagnosis	Cat.No.
AIV Detection Kit	AIV type A (Matrix gene)	PD55-01
AIV H5 Detection Kit	Subtype H5	PD55-02
AIV H7 Detection Kit	Subtype H7	PD55-03
AIV H9 Detection Kit	Subtype H9	PD55-04

PCR Q & A

Why do AIV PCR Kits require for AIV diagnosis?

- Because the culture methods are time consuming and heavily contaminated (false negative)
- Because the serological methods are not specific enough (false positive) or when the delay of seroconversion is long

What are the samples for PCR assays?

- Sample: RNA extracted from blood, tissue, secretion of birds, horse, pig or human

How accurate?

- Sensitive: 5 fg of viral RNA detectable, detects early at onset of disease as early as 1~2 days
- Precise: Highly specific primers detect only avian influenza viruses.

Advantage

- More sensitivity than general PCR
- Simple and fast PCR procedures

Avian General PCR Kits NDV, IBDV detection Kit



Background

BIONOTE NDV and IBDV Detection Test Kit is a complete, easy PCR kit for detection of Specific protein gene from Newcastle disease virus and Infectious Bursal Disease virus. These detection Kit are very accurate and high efficient, and RT-reaction can be done from 50fg to 500ng of template RNA. It is developed with the best condition of synthesis first-strand cDNA, so it is useful for check a low copy of DNA transcription.

Specifications

- ✓ High Sensitivity and Specificity : BIONOTE has specific Primer and Probes system
- ✓ Excellent reproducibility
- ✓ Easy to use: BioNote Detection Kits contain all components for the detection of NDV or IBDV.
- ✓ Saving labor costs& time in the laboratory
- ✓ Includes positive and negative control

NDV or IBDV Detection Kits Procedure




1. Sampling



2. RNA Extraction



3. Amplification



4. Detection

Materials in NDV or IBDV Detection kit (96tests/Kit)

Contents	Volume
Negative Control	0.5ml
Positive Control	0.5ml
20X RT Mix	0.1ml
20X PCR Enzyme Mix	0.1ml
Detection Solution	0.8ml

Type of NDV, and IBDV Real-Time detection Kit

Product	Diagnosis	Cat.No.
Lento NDV Detection Kit	Lentogenic Newcastle Disease virus	PD55-07
Velo NDV Detection Kit	Velogenic Newcastle Disease virus	PD55-08
IBDV Detection Kit	Infectious Bursal Disease virus	PD55-11

Avian Real-Time PCR Kit AIV detection Kit



PCR

Background

AIV Real-Time Detection Kit from BIONOTE is useful for detecting the infection of avian influenza A viruses. The kit can be exactly performed to detect RNA of avian influenza A virus, H5, H7, H9 according to subtype, so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary type.

Specifications

- ✓ **High Sensitivity** : Each BIONOTE Real-Time Detection Kit provides highly effective pathogen detection in a sample material.
→ Possible detection up to 10 copies/ul
- ✓ **High Specificity**: BIONOTE owing optimal & unique primers and probes system guarantees a highly specific result
- ✓ **Excellent reproducibility**: One step Real-Time PCR Kit using Taqman probe
- ✓ **Quantitative and qualitative analysis.**
- ✓ **Complete result within 2 hours** for after RNA extraction
- ✓ **Higher detection limit than general PCR** : 100~1,000 times more sensitive
- ✓ **Simple test procedure**: Adding only detection solution and enzyme
- ✓ **BioNote AIV Detection Kits for AIV** contain all component for detection from sample.
- ✓ **Feces, cloaca, and organ homogenates samples** are recommended

Materials in AIV Detection Kit (96tests/Kit)

Contents	Volume
Standard 1 (1x10 ⁴ copies/μl)	200μl x 1
Standard 2 (1x10 ⁵ copies/μl)	200μl x 1
Standard 3 (1x10 ⁶ copies/μl)	200μl x 1
Standard 4 (1x10 ⁷ copies/μl)	200μl x 1
Detection Solution*1	460μl x 1
Nuclease Free buffer	1500μl x 1
2X Enzyme buffer	1000μl x 1
Enzyme Mix	40μl x 1
Rox Reference Dye *2	40μl x 1

*1: Probe is labeled at the 5'-end with the reporter molecule 6-carboxyfluorescein (FAM).

*2: The kits is provided extra Rox dye for ABI real-time PCR

Real-Time Detection Kits Procedure

Extraction the template RNA

Prepare the PCR MiX

React the testing sample, Standard 1,2,3,4 and negative control at each test tube

Set up the cycle program as following the Applicable Real-Time PCR machine

Set up the report dye as FAM

List of AIV Real-Time Detection Kit

Product	Diagnosis	Cat.No.
AIV A Real -Time Detection Kit	All type A (Matrix gene)	PD65-01
AIV H5 Real -Time Detection Kit	Subtype H5	PD65-02
AIV H7 Real -Time Detection Kit	Subtype H7	PD65-03
AIV H9 Real -Time Detection Kit	Subtype H9	PD65-04
AIV H5/N1 Real-Time Detection Kit	Subtype H5 N1	PD65-016
AIV N1 Real-Time Detection Kit	Subtype N1	PD65-017
AIV H5/N1 Duplex Detection Kit	Subtype H5 and N1	PD65-051

Avian Real-Time PCR Kit Other poultry disease virus detection Kit



Background

Other poultry disease virus Real-Time Detection Kit from BIONOTE is useful for detecting the infection of REV, MDV, NDV, IBV, TRT, IBDV, REO, ORT, Myco or ALV. The kit can be exactly performed to detect RNA or DNA of each virus according to disease, so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary type.

PCR

Specifications

- ✓ **High Sensitivity** : Each BIONOTE Real-Time Detection Kit provides highly effective pathogen detection in a sample material.
→ Possible detection up to 10 copies/ul
- ✓ **High Specificity**: BIONOTE owing optimal & unique primers and probes system guarantees a highly specific result
- ✓ **Excellent reproducibility**: One step Real-Time PCR Kit using Taqman probe
- ✓ **Quantitative and qualitative analysis.**
- ✓ **Complete result within 2 hours** for after RNA or DNA extraction
- ✓ **Higher detection limit than general PCR** : 100~1,000 times more sensitive
- ✓ **Simple test procedure**: Adding only detection solution and enzyme
- ✓ **BioNote Detection Kits for each virus** contain all component for detection from sample.

Specifications

- Reaction Buffer
- dNTP
- Sense Primer
- Anti Sense Primer
- Taqman probe
- UNG
- Rox reference dye

- Easy handling
- Fast procedure
- Low contamination risk
- High reproducibility

List of Other poultry disease Real-Time Detection Kit

Product	Diagnosis	Cat.No.
REV Real-Time Detection Kit	Reticuloendotheliosis Virus	PD65-05
MDV Real-Time Detection Kit	Marek's disease virus	PD65-06
Lento NDV Real-Time Detection Kit	low pathogenic NDV virus (Lentogenic)	PD65-07
Velo NDV Real-Time Detection Kit	high pathogenic NDV virus (Velogenic)	PD65-08
IBV Real-Time Detection Kit	Infectious Bronnchitis Virus	PD65-09
TRT Real-Time Detection Kit	Infectious Turkey Rhinotracheitis	PD65-10
IBDV Real-Time Detection Kit	Infectious Bursal Disease Virus	PD65-11
REO Real-Time Detection Kit	Avian Reovirus	PD65-12
ORT Real-Time Detection Kit	Ornithobacterium rhinotracheale	PD65-13
MYCO G/S Real-Time Detection Kit	Avian Mycoplasma gallisepticum, and M. iowae, M. meleagridis , M. synoviae	PD65-14
ALV Real-Time Detection Kit	Avian Leukosis Virus	PD65-15

Swine

Real-Time PCR Kit

PRRSV, PCV2, CSFV detection Kit



PCR

Background

***PCV-2 Real -Time Detection Kit and *CSFV Real -Time Detection Kit from BIONOTE are** useful for detection the infection of Porcine Circovirus or Classical swine fever virus. The Kit can be exactly performed to detect DNA of Porcine Circovirus or Classical swine fever virus, so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary tubes

***PRRSV Real -Time Detection Kit from BIONOTE is** useful for **differential** detection the infection of EU/US strain of Porcine Reproductive and Respiratory Syndrome Virus. The Kit can be exactly performed to detect DNA of Porcine Reproductive and Respiratory Syndrome Virus (EU/US), so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary tubes


PCV-2 or CSFV Real-Time Detection Kit

Contents	Volume
Standard 1 (1x10 ⁴ copies/μl)	200μl x1
Standard 2 (1x10 ⁵ copies/μl)	200μl x1
Standard 3 (1x10 ⁶ copies/μl)	200μl x1
Standard 4 (1x10 ⁷ copies/μl)	200μl x1
Detection Solution PCV-2	500μl x1
Nuclease Free Water	1500μl x1
2X Master Mix	1000μl x1
Rox Reference Dye	40μl x1

PRRSV Real-Time Detection Kit

Contents	Volume
Standard 1 (1x10 ⁴ copies/μl)	200μl x1
Standard 2 (1x10 ⁵ copies/μl)	200μl x1
Standard 3 (1x10 ⁶ copies/μl)	200μl x1
Standard 4 (1x10 ⁷ copies/μl)	200μl x1
Detection Solution PRRS EU	460μl x1
Detection Solution PRRS US	460μl x1
Nuclease Free Water	1500μl x1
2X Enzyme buffer	1000μl x1
Enzyme Mix	40μl x1
Rox Reference Dye	40μl x1

Swine Real-Time Detection Kits Procedure



List of PCV and PRRSV Real-Time Detection Kit

Product	Diagnosis	Cat.No.
PCV-2 Real-Time Detection Kit	Porcine Circovirus	PD64-03
PRRSV Real-Time Detection Kit	Porcine Reproductive and Respiratory Syndrome Virus (Differential diagnosis of EU/US strain)	PD64-04
CSFV Real-Time Detection Kit	Classical swine fever virus	PD64-05

Swine Real-Time PCR Kit SIV Detection Kit



Background

*SIV Real -Time Detection Kit from BIONOTE is useful for detection the infection of Swine Influenza Virus. There are 2 types of SIV Real-Time Detection Kit. One is Swine Influenza Virus (Influenza A H1N1) detection Kit and the other distinguishes Classical Swine Influenza virus between and Novel Swine Influenza virus. These Kits can be exactly performed to detect RNA of Swine Influenza virus according to subtype, so it can be used for both qualitative and quantitative analysis. It works most of Real-Time PCR apparatuses of block and capillary tubes.

PCR

Specifications

- High Sensitivity : Each BIONOTE Real-Time Detection Kit provides highly effective pathogen detection in a sample material.
→ Possible detection up to 10 copies/ul
- High Specificity: BIONOTE owing optimal& unique primers and probes system guarantees a highly specific result
- Differential Diagnosis: Original Swine Influenza and Novel Influenza (H1N1) are easily distinguished by SIV Real-Time Detection Kit (PD64-02)
- Excellent reproducibility: One step Real-Time PCR Kit using Taqman probe
- Quantitative and qualitative analysis.
- Complete result within 2 hours for after RNA extraction
- Higher detection limit than general PCR : 100~1,000 times more sensitive

SIV Real -Time Detection Kit

Contents	Volume
Standard 1 (1x10 ⁴ copies/μl)	200μlx1
Standard 2 (1x10 ⁵ copies/μl)	200μlx1
Standard 3 (1x10 ⁶ copies/μl)	200μlx1
Standard 4 (1x10 ⁷ copies/μl)	200μlx1
Detection Solution Swine Influenza virus	500μlx1
Nuclease Free Water	1500μlx1
2X Master Mix	1000μlx1
Rox Reference Dye	40μlx1

Influenze A(H1N1)-Swine Real-Time Detection Kit

Contents	Volume
Positive control	200μlx1
Detection Solution Influenza A (M gene)	460μlx1
Detection Solution Novel Influenza A (H1 gene)	460μlx1
Nuclease Free Water	1500μlx1
2X Enzme buffer	1000μlx1
Enzyme Mix	80μlx1
Rox Reference Dye	800μlx1

List of SIV Real-Time detection Kit






Product	Diagnosis	Cat.No.
Influenza A(H1N1)-Swine Real-Time Detection Kit	Novel Influenza A (H1n1)	PD64-01
SIV Real-Time Detection Kit	Original Swine Influenza virus	PD64-02

Real-Time Fluorescence Detection System





CURRENT AVAILABLE PRODUCT LIST


Rapid

	Species	Product	Cat. No.	Type	Packing Size	Description	Sample
PET Animal		Rapid CPV Ag Test Kit	RG11-01	Device	1Test x 10/Kit	Detection of Canine Parvovirus antigen	Feces
		Rapid Canine Heartworm Ag Test Kit	RG11-02	Device	1Test x 10/Kit	Detection of Canine Heartworm antigen	Whole blood, plasma or serum
		Rapid CDV Ag Test Kit	RG11-03	Device	1Test x 10/Kit	Detection of Canine Distemper virus antigen	Conjunctiva, feces, saliva, urine, plasma, or serum
		Rapid CCV Ag Test Kit	RG11-04	Device	1Test x 10/Kit	Detection of Canine Coronavirus antigen	Feces
		Rapid CPV Ag/CCV Ag Test Kit	RG11-05	Dual device	1Test x 10/Kit	Detection of Parvovirus antigen/Coronavirus antigen	Feces
		Rapid CIV Ag Test Kit	RG11-07	Device	1Test x 10/Kit	Detection of Canine influenza virus antigen	Conjunctiva and nasal swab
		Rapid CDV Ag/CAV Ag Test Kit	RG11-08	Dual device	1Test x 10/Kit	Detection of Distemper antigen/Canine adenovirus antigen	Conjunctiva and nasal swab
		Rapid CIRCD-3 Ag Test Kit	RG11-09	Dual device	1Test x 10/Kit	Detection of Distemper antigen/Canine adenovirus/Canine influenza virus antigen	Conjunctiva and nasal swab
		Rapid CDV Ag/CIV Ag Test Kit	RG11-10	Dual device	1Test x 10/Kit	Detection of Distemper antigen/Canine Influenza antigen	Conjunctiva and nasal swab
		Rapid CPV/CCV/Giardia Ag Test Kit	RG11-12	Triple Device	1Test x 5/Kit	Detection of Parvovirus antigen/Coronavirus antigen/ Giardia antigen	Feces
		Rapid Rabies Ag Test Kit	RG18-01	Device	1Test x 10/Kit	Detection of Rabies virus antigen	Saliva & brain homogenates
		Rapid Rota Ag Test Kit	RG18-03	Device	1Test x 10/Kit	Detection of Rotavirus antigen in bovine, porcine, canine	Diarrheal faeces
		Rapid Giardia Ag Test kit	RG18-04	Device	1Test x 10/Kit	Detection of Giardia antigen	Feces
		Rapid CPV Ab Test Kit	RB21-01	Device	1Test x 10/Kit	Titration of Parvovirus antibody for vaccination & prognosis	Whole blood, plasma or serum
	Rapid CDV Ab Test Kit	RB21-02	Device	1Test x 10/Kit	Titration of Distemper virus antibody for vaccination & prognosis	Blood	
	Rapid Canine Brucella Ab Test Kit	RB21-03	Device	1Test x 10/Kit	Detection of Brucella canis antibody	Whole blood, plasma or serum	
	Rapid Leishmania Ab Test Kit	RB21-04	Device	1Test x 10/Kit	Detection of Leishmania antibody	Whole blood, plasma or serum	
		Rapid FeLV Ag Test Kit	RG12-01	Device	1Test x 10/Kit	Detection of Feline Leukemia virus antigen	Whole blood, plasma or serum
		Rapid FPV Ag Test Kit	RG12-03	Device	1Test x 10/Kit	Detection of Feline Panleukopenia virus antigen	Feces
		Rapid FIV Ab/FeLV Ag Test Kit	RC12-04	Dual Device	1Test x 10/Kit	Detection of Feline Leukemia virus antigen and Feline Immunodeficiency virus antibody	Whole blood, plasma or serum
Rapid FCoV Ag Test Kit		RG12-05	Device	1Test x 10/Kit	Detection of Feline Corona virus antigen	Feces	
Rapid FIV Ab Test Kit		RB22-01	Device	1Test x 10/Kit	Detection of Feline Immunodeficiency virus antibody	Whole blood, plasma or serum	
Rapid FCoV Ab Test Kit		RB22-03	Device	1Test x 10/Kit	Detection of Feline Corona virus antibody	Whole blood, plasma or serum	
Industrial Animal		Rapid Rota Ag Test Kit	RG18-03	Device	1Test x 10/Kit	Detection of Rotavirus antigen in bovine, porcine, canine	Diarrheal faeces
		Rapid Bovine Brucella Ab Test Kit	RB23-01	Device	1Test x 10/Kit	Detection of Brucella abortus antibody	Whole blood, plasma, serum or milk
		Rapid Bovine TB Ab Test Kit	RB23-02	Device	1Test x 10/Kit	Detection of Mycobacterium bovis antibody	Plasma, serum
		Rapid FMD NSP Ab Test Kit	RB28-02	Device	1Test x 10/Kit	Detection of Foot and mouth disease virus antibody field infected	Whole blood, plasma or serum
		Rapid PED Ag Test Kit	RG14-01	Device	1Test x 10/Kit	Detection of Porcine Epidemic Diarrhea virus antigen	Diarrheal faeces
		Rapid TGE Ag Test Kit	RG14-02	Device	1Test x 10/Kit	Detection of Transmissible Gastroenteritis virus antigen	Diarrheal faeces
		Rapid TGE Ag/PED Ag Test Kit	RG14-03	Dual device	1Test x 10/Kit	Detection of TGE Ag and PED Ag	Diarrheal faeces
		Rapid PED Ag/Rota Ag Test Kit	RG14-04	Dual device	1Test x 10/Kit	Detection of PED antigen and Rotavirus antigen in pig	Diarrheal faeces
		Rapid Rota Ag Test Kit	RG18-03	Device	1Test x 10/Kit	Detection of Rotavirus antigen in bovine, porcine, canine	Diarrheal faeces
		Rapid FMD NSP Ab Test Kit	RB28-02	Device	1Test x 10/Kit	Detection of Foot and mouth disease virus antibody field infected	Whole blood, plasma or serum
		Rapid SIV Ag Test Kit	RG14-05	Device	1Test x 10/Kit	Detection of Swine Influenza virus antigen	Nasal fluid or tracheal swab and lung tissue
		Rapid AIV Ag Test Kit	RG15-01	Multidevice	10Tests x 3/Kit	Detection of Avian influenza type A virus antigen	Diarrhoea faeces, Cloaca Feces, or Trachea swab
		Rapid NDV Ag Test Kit	RG15-03	Device	1Test x 10/Kit	Detection of Newcastle Disease virus antigen	Cloaca Feces, Trachea
		Rapid IBDV Ag Test Kit	RG15-04	Device	1Test x 10/Kit	Detection of Infectious Bursal Disease virus antigen	Bursa of Fabricious, Feces
		Rapid H5 AIV Ag Test Kit	RG15-05	Device	1Test x 25/Kit	Detection of avian influenza type A subtype H5 virus antigen	Diarrhoea faeces, Cloaca Feces, or Trachea swab
Rapid AIV type A/H5 AIV Ag Test Kit		RG15-06	Dual device	1Test x 10/Kit	Detection of avian influenza type A and H5 antigen	Diarrhoea faeces, Cloaca Feces, or Trachea swab	
Rapid AIV Ag/NDV Ag Test Kit		RG15-09	Dual device	1Test x 10/Kit	Detection of avian influenza type A and NDV antigen	Diarrhoea faeces, Cloaca Feces, or Trachea swab	
Rapid AIV Ab Test Kit		RB25-01	Multidevice	10Tests x 3/Kit	Detection of Avian Influenza Virus Type A Antibody	Plasma, serum	



ELISA

	Species	Product	Cat. No.	Type	Packing Size	Description	Sample
Industrial Animal		Bovine Brucella Ab ELISA	EB43-01	Microplate	480 wells/Kit	Detection of Brucella abortus antibody	Serum, Plasma, Raw milk
			EB43-02		960 wells/Kit		
		BTB Ab ELISA	EB43-03	Microplate	480 wells/Kit	Detection of Mycobacterium bovis antibody	Serum
		FMD NSP Ab ELISA	EB48-01	Microplate	480 wells/Kit	Detection of Foot-and-mouth disease virus antibody field infected in cattle, pig, goat, sheep	Serum, Plasma
		CSFV Ab ELISA	EB44-13	Microplate	480 wells/Kit	Detection of Classical swine fever virus antibody for routine screening	Serum, Plasma
	PRRS Ab ELISA	EB44-04	Microplate	480 wells/Kit	Detection of Porcine Reproductive and Respiratory Syndrome (PRRS) antibody	Serum, Plasma	




ELISA

	Species	Product	Cat. No.	Type	Packing Size	Description	Sample
Industrial Animal		AIV Ab ELISA	EB45-02	Microplate	480 wells/Kit	Detection of AIV type A antibody in chicken, duck, goose.	Serum
			EB45-52		96 wells/Kit		
		H5 AIV Ab ELISA	EB45-03	Microplate	480 wells/Kit	Detection of AIV subtype H5 antibody in chicken, duck, goose.	Serum
			EB45-53		96 wells/Kit		
		NDV Velo Ab ELISA	EB45-01	Microplate	480 wells/Kit	Detection of NDV velogenic virus antibody in chicken	Serum
			EB45-51		96 wells/Kit		

Conventional PCR

	Species	Product	Cat. No.	Type	Packing Size	Description	Sample
Industrial Animal		AIV Detection Test Kit	PD55-01	Conventional PCR	96Tests/kit	Detection of AIV type A	sercretion, organ, blood, feces, fluid etc.
		H5 AIV Detection Test Kit	PD55-02	Conventional PCR	96Tests/kit	Detection of AIV subtype H5	sercretion, organ, blood, feces, fluid etc.
		H7 AIV Detection Test Kit	PD55-03	Conventional PCR	96Tests/kit	Detection of AIV subtype H7	sercretion, organ, blood, feces, fluid etc.
		H9 AIV Detection Test Kit	PD55-04	Conventional PCR	96Tests/kit	Detection of AIV subtype H9	sercretion, organ, blood, feces, fluid etc.
		Lento NDV Detection Test Kit	PD55-07	Conventional PCR	96Tests/kit	Detection of low pathogenic NDV virus (Lentogenic)	sercretion, organ, blood, feces, fluid etc.
		Velo NDV Detection Test Kit	PD55-08	Conventional PCR	96Tests/kit	Detection of high pathogenic NDV virus (Velogenic)	sercretion, organ, blood, feces, fluid etc.
		IBDV Detection Test Kit	PD55-11	Conventional PCR	96Tests/kit	Detection of Infectious Bursal Disease Virus	sercretion, organ, blood, feces, fluid etc.
PET Animal		CDV Detection Test Kit	PD 51-03	Conventional PCR	25Tests/kit	Detection of Canine Desteemper Virus	Conjunctiva, feces, saliva, urine, plasma, or serum
		Leishmaniasis Detection Test Kit	PD 51-04	Conventional PCR	25Tests/kit	Detection of Leishmania	Whole blood, plasma or serum
		E.canis Detection Test Kit	PD 51-05	Conventional PCR	25Tests/kit	Detection of E.canis	Whole blood, plasma or serum

Real-Time PCR

	Species	Product	Cat. No.	Type	Packing Size	Description	Sample
Industrial Animal		AIV A Real-Time Detection Test Kit	PD65-01	Real-Time PCR	96Tests/kit	Detection of AIV type A	sercretion, organ, blood, feces, fluid etc.
		H5 AIV Real-Time Detection Test Kit	PD65-02	Real-Time PCR	96Tests/kit	Detection of AIV subtype H5	sercretion, organ, blood, feces, fluid etc.
		H7 AIV Real-Time Detection Test Kit	PD65-03	Real-Time PCR	96Tests/kit	Detection of AIV subtype H7	sercretion, organ, blood, feces, fluid etc.
		H9 AIV Real-Time Detection Test Kit	PD65-04	Real-Time PCR	96Tests/kit	Detection of AIV subtype H9	sercretion, organ, blood, feces, fluid etc.
		REV Real-Time Detection Test Kit	PD65-05	Real-Time PCR	96Tests/kit	Detection of Reticuloendotheliosis Virus	sercretion, organ, blood, feces, fluid etc.
		MDV Real-Time Detection Test Kit	PD65-06	Real-Time PCR	96Tests/kit	Detection of Marek's disease virus	sercretion, organ, blood, feces, fluid etc.
		Lento NDV Real-Time Detection Test Kit	PD65-07	Real-Time PCR	96Tests/kit	Detection of low pathogenic NDV virus (Lentogenic)	sercretion, organ, blood, feces, fluid etc.
		Velo NDV Real-Time Detection Test Kit	PD65-08	Real-Time PCR	96Tests/kit	Detection of high pathogenic NDV virus (Velogenic)	sercretion, organ, blood, feces, fluid etc.
		IBV Real-Time Detection Test Kit	PD65-09	Real-Time PCR	96Tests/kit	Detection of Infectious Bronchitis Virus	sercretion, organ, blood, feces, fluid etc.
		TRT Real-Time Detection Test Kit	PD65-10	Real-Time PCR	96Tests/kit	Detection of Infectious Turkey Rhinotracheitis	sercretion, organ, blood, feces, fluid etc.
		IBDV Real-Time Detection Test Kit	PD65-11	Real-Time PCR	96Tests/kit	Detection of Infectious Bursal Disease Virus	sercretion, organ, blood, feces, fluid etc.
		REO Real-Time Detection Test Kit	PD65-12	Real-Time PCR	96Tests/kit	Detection of Avian Reovirus	sercretion, organ, blood, feces, fluid etc.
		ORT Real-Time Detection Test Kit	PD65-13	Real-Time PCR	96Tests/kit	Detection of Orthobacterium rhinotracheale	sercretion, organ, blood, feces, fluid etc.
		MYCO G/S Real-Time Detection Test Kit	PD65-14	Real-Time PCR	96Tests/kit	Detection of Avian Mycoplasma gallisepticum, and M. iowae, M. meleagridis, M. synoviae	sercretion, organ, blood, feces, fluid etc.
		ALV Real-Time Detection Test Kit	PD65-15	Real-Time PCR	96Tests/kit	Detection of Avian Leukosis Virus	sercretion, organ, blood, feces, fluid etc.
		AIV H5/N1 Real-Time Detection Test Kit	PD65-16	Real-Time PCR	96Tests/kit	Detection of AIV H5N1 (Only FAM)	sercretion, organ, blood, feces, fluid etc.
		AIV N1 Real-Time Detection Test Kit	PD65-17	Real-Time PCR	96Tests/kit	Detection of AIV N1	sercretion, organ, blood, feces, fluid etc.
	AIV H5/N1 Duplex Real-Time Detection Test Kit	PD65-51	Real-Time PCR	96Tests/kit	Detection of AIV H5 and N1	sercretion, organ, blood, feces, fluid etc.	
		Influenza A(H1N1)-Swine Real-Time Detection Test Kit	PD64-01	Real-Time PCR	96Tests/kit	Detection of Novel Influenza A(H1N1) Virus (Swine)	sercretion, organ, blood, feces, fluid etc.
		SIV Real-Time Detection Test Kit	PD64-02	Real-Time PCR	96Tests/kit	Detection of Swine Influenza Virus	sercretion, organ, blood, feces, fluid etc.
PCV-2 Real-Time Detection Test Kit		PD64-03	Real-Time PCR	96Tests/kit	Detection of Porcine Circovirus	sercretion, organ, blood, feces, fluid etc.	
PRRSV Real-Time Detection Test Kit		PD64-04	Real-Time PCR	96Tests/kit	Detection of Porcine Reproductive and Respiratory syndrome Virus	sercretion, organ, blood, feces, fluid etc.	
CSFV Real-Time Detection Test Kit		PD64-05	Real-Time PCR	96Tests/kit	Detection of Classical swine fever virus	sercretion, organ, blood, feces, fluid etc.	
PET Animal		CHW Real-Time Detection Kit	PD61-02	Real-Time PCR	96Tests/kit	Detection of Canine Heartworm	
		Influenza A(H1N1) Real-Time Detection Test Kit (2 Wells type)	PD65-18	Real-Time PCR	96Tests/kit	Detection of Novel Influenza A(H1N1) Virus (Human)	sercretion, organ, blood, fluid etc.
		Influenza A(H1N1) Real-Time Detection Test Kit (4 Wells type)	PD65-19	Real-Time PCR	96Tests/kit	Detection of Novel Influenza A(H1N1) Virus (Human)	sercretion, organ, blood, fluid etc.
		OneStep RT-PCR Kit	RK50-01	RT-PCR	100 Reactions		sercretion, organ, blood, fluid etc.
		SYBR Green RT-PCR Kit	RK50-02	RT-PCR	1000 Reactions		sercretion, organ, blood, fluid etc.
		Viral RNA purification Kit	RK50-90		250 Reactions		sercretion, organ, blood, fluid etc.





BIONOTE Product catalog

Third Edition CIA03-03

Address

2-9 Seogu-dong, Hwaseong-si,
Gyeonggi-do, Korea 445-170

Tel. 82 31 211 0516
Fax. 82 31 8003 0618

www.bionote.co.kr